



Touch Prime Ultrasound System

DICOM Conformance Statement

2017-05-18

Document # 9J9407

The information contained herein is based on the experience and knowledge relating to the subject matter gained by Carestream Health, Inc. prior to publication.

No patent license is granted by this information.

Carestream Health reserves the right to change this information without notice, and makes no warranty, express or implied, with respect to this information. Carestream Health shall not be liable for any loss or damage, including consequential or special damages, resulting from any use of this information, even if loss or damage is caused by Carestream Health's negligence or other fault.

Network Services

Networking SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Transfer		
Ultrasound Image Storage	Yes	No
Ultrasound Multiframe Image Storage	Yes	No
Storage Commitment Push Model	Yes	No
Comprehensive SR	Yes	No
Query/Retrieve		
Patient Root Q/R - FIND	Yes	No
Patient Root/Study Root Q/R - MOVE	Yes	No
Workflow Management		
Modality Worklist	Yes	No
Modality Performed Procedure Step	Yes	No
Print Management		
Basic Grayscale Print Management	Yes	No
Basic Color Print Management	Yes	No

Supported Structured Report Templates

Concept Name	Template ID
Vascular Ultrasound Procedure Report	TID 5100
Structured Report	TID CS1000

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Contents

1	Introduction.....	7
1.1	Scope and Field of Application	7
1.2	Intended Audience	8
1.3	Used Definitions, Terms, Symbols and Abbreviations	8
1.4	References	10
1.5	Important Note to the Reader.....	10
	1.5.1 Interoperability.....	10
	1.5.2 Validation	10
	1.5.3 New Versions of the DICOM Standard.....	11
2	Implementation Model - Network	12
2.1	Application Data Flow Diagram	12
	2.1.1 Verify Communication.....	14
	2.1.2 Store Image and Reports.....	14
	2.1.3 Storage Commitment	14
	2.1.4 Retrieve HIS/RIS Worklist	15
	2.1.5 Performed Procedure Step	15
	2.1.6 Print Image	15
	2.1.7 Query/Retrieve.....	15
2.2	Functional Definition of Application Entity	15
	2.2.1 Real World Activity: Verification	16
	2.2.2 Real World Activity: Storage	16
	2.2.3 Real World Activity: Request Worklist	16
	2.2.4 Real World Activity: Performed Procedure Step	17
	2.2.5 Real World Activity: Print	17
	2.2.6 Real World Activity: Query/Retrieve.....	17

2.3	Sequencing of Real-World Activities	18
3	AE Specifications.....	19
3.1	AE Specification	19
3.1.1	Association Establishment Policies	19
3.1.2	Association Initiation Policy	20
3.1.3	Association Acceptance Policy	36
4	Communication Profiles	37
4.1	TCP/IP Stack.....	37
4.1.1	Physical Media Support.....	37
5	Extensions / Specializations / Privatizations	38
6	Configuration	39
6.1	AE Title/Presentation Address Mapping	39
6.2	Configurable Parameters.....	39
6.2.1	Local AE	39
6.2.2	Remote AE.....	39
6.2.3	Storage	39
6.2.4	Storage Commitment	40
6.2.5	Printing	40
6.2.6	Query/Retrieve.....	40
7	Support of Extended Character Sets	41
8	Security Profiles	41
9	Appendix: Overview of the Applied Information Object Definitions	43
9.1	Overview of the Applied Ultrasound (US) Storage SOP Classes.....	43
9.2	Overview of the Applied Modality Worklist IOD.....	50
9.3	Overview of the Applied Basic Grayscale Print Management Meta IOD.....	53
9.3.1	Overview of the Applied Basic Film Session IOD.....	53
9.3.2	Overview of the Applied Basic Film Box IOD.....	53

9.3.3	Overview of the Applied Basic Grayscale Image Box IOD.....	54
9.4	Overview of the Applied Basic Color Print Management Meta IOD.....	54
9.4.1	Overview of the Applied Basic Film Session IOD.....	54
9.4.2	Overview of the Applied Basic Film Box IOD.....	55
9.4.3	Overview of the Applied Basic Color Image Box IOD.....	56
9.5	Overview of the Applied Modality Performed Procedure Step IOD.....	56
9.6	Overview of the Applied Storage Commitment IOD.....	59
9.7	Overview of the Applied Comprehensive SR IOD.....	60
9.7.1	Overview of SR Document Content Descriptions.....	63
9.8	Overview of the Applied Patient Root Query IOD.....	74
9.9	Overview of the Applied Patient Root/Study Root Retrieve IOD.....	105

1 Introduction

This section provides general information about the purpose, scope and contents of this Conformance Statement.

This document describes the conformance to the DICOM 3.0 Standard of the Touch Prime Ultrasound System.

The Touch Prime Ultrasound System is a device that:

- Generates ultrasound images and other data that:
 - Can be sent using DICOM standard protocols and definitions to network archive servers.
 - Can be printed on a remote printer.
- Generates ultrasound comprehensive structure reporting and other data that:
 - Can be sent using DICOM standard protocols and definitions to network archive servers.
- Can retrieve data from a Radiology Information System (RIS).
- Can notify the Remote Modality Performed Procedure Step server on the procedures performed.
- Can query/retrieve ultrasound images and structured reports with a set of key attributes.

1.1 Scope and Field of Application

The scope of this DICOM Conformance Statement is to facilitate data exchange with equipment from Carestream . This document specifies the compliance to the DICOM standard, formally called the DICOM PS 3.X standards. It contains a short description of the applications involved and provides technical information about the data exchange capabilities of the equipment. The main elements describing these capabilities are: the supported DICOM Service Object Pair (SOP) Classes, Application Profiles, Roles, Information Object Definitions (IOD) and Transfer Syntaxes.

The field of application is the integration of the Carestream equipment into an environment of medical devices.

This Conformance Statement should be read in conjunction with the DICOM standard and its addenda.

1.2 Intended Audience

This Conformance Statement is intended for:

- (Potential) customers.
- System integrators of medical equipment.
- Marketing staff interested in system functionality.
- Software designers implementing DICOM interfaces.

It is assumed that the reader is familiar with the DICOM standard.

1.3 Used Definitions, Terms, Symbols and Abbreviations

DICOM definitions, terms and abbreviations are used throughout this Conformance Statement. For a description of these, see DICOM PS 3.3 and PS 3.4.

The following symbols and abbreviations are used in this document:

AE	Application Entity
AP	Application Profile
BCID	Baseline Context Group Identifier
BKCMR	Content Mapping Resource
DCID	Defined Context Group Identifier
DCMR	DICOM Content Mapping Resource
DICOM	Digital Imaging and Communications in Medicine
DIMSE	DICOM Message Service Element
DIMSE-C	DICOM Message Service Element-Composite
DIMSE-N	DICOM Message Service Element-Normalized
EV	Enumerated Value
HIS	Hospital Information System

IE	Information Entity
IFF	If and only if
IOD	Information Object Definition
ISO	International Standards Organization
MPPS	Modality Performed Procedure Step
NEMA	National Electrical Manufacturers Association
PACS	Picture Archiving and Communication System
PDU	Protocol Data Unit
Q/R	Query/Retrieve
RIS	Radiology Information System
RWA	Real-World Activity
SC	Service Class
SCP	Service Class Provider
SCU	Service Class User
SNOMED	Systematized Nomenclature of Medicine
SOP	Service-Object Pair
SR	Structured Reporting
TCP/IP	Transmission Control Protocol/Internet Protocol
TID	Template ID
UCUM	Unified Code for Units of Measure
UID	Unique Identifier
US	Ultrasound
VR	Value Representation

The following upper-case abbreviations represent specific Attributes:

CM	Code Meaning (0008, 0104)
CSD	Coding Scheme Designator (0008, 0102)
CSV	Coding Scheme Version (0008,0103)
CV	Code Value (0008, 0100)

1.4 References

DICOM The Digital Imaging and Communications in Medicine (DICOM) standard:
DICOM PS 3.X (X refers to the part 1 - 20 and Supplements).
National Electrical Manufacturers Association (NEMA)
1300 N. 17th Street, Suite 900, Rosslyn, Virginia 22209, United States of America.
Can also be downloaded from <http://medical.nema.org/>

1.5 Important Note to the Reader

This Conformance Statement by itself does not guarantee successful interoperability of Carestream equipment with non-Carestream equipment. The user should be aware of the following issues:

- Interoperability.
- Validation.
- New (or old) versions of the DICOM standard.

1.5.1 Interoperability

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of Carestream equipment with non-Carestream equipment.

It is the user's responsibility to thoroughly analyze the application requirements and to specify a solution that integrates Carestream equipment with non-Carestream equipment.

1.5.2 Validation

Carestream equipment has been tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.

Where Carestream equipment is linked to non-Carestream equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation tests will be necessary to ensure the functionality, performance, accuracy and stability of image, report and image or report related data.

It is the user's responsibility to specify the appropriate test and to carry out the additional validation tests.

1.5.3 New Versions of the DICOM Standard

The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. Carestream plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, Carestream reserves the right to make changes to its products or to discontinue their delivery.

2 Implementation Model - Network

Carestream's Ultrasound Systems with DICOM option activated has implemented the DICOM functionality as one Application Entity.

The AE contains the following DICOM functionality:

- Verify communication to a remote AE.
- Retrieve a worklist for ultrasound modality from a Radiology Information System (RIS) or a Hospital Information System (HIS).
- Send Step by Step information on the procedures performed to a remote server.
- Transfer ultrasound images to a remote storage system.
- Transfer ultrasound comprehensive structured reports to a remote storage system.
- Transfer ultrasound video clips as multi-frame images to a remote storage system.
- Request Storage Commitment of images and reports to the remote storage system.
- Print ultrasound images on a remote printer.
- Query information about ultrasound images and structured reports from a remote storage system.
- Retrieve ultrasound images and structured reports from a remote storage system.

2.1 Application Data Flow Diagram

The Implementation Model for the AE, the DICOM networking service for the Touch Prime Ultrasound System, is depicted in Figure 2-1.




Figure 2-1 AE Implementation Model

The AE can send image and report storage objects. It receives requests from the operator to transmit an image and report to a specific DICOM destination. It receives the Storage Commitment reports from the remote storage server and passes it to the Modality. It can query for information from external sources. Worklist demographic queries can be initiated and executed by the AE. It can query for patient, study, series, and composite object instance data. It can issue instance retrieval requests and receive the retrieved instances. It receives the performed procedure step notifications from the Modality and initiates an association and transmits the information to the remote MPPS Server. The operator can

initiate an association with a DICOM-compliant grayscale or color printer to print a selected image.

The AE is implemented as a Windows®-based application.

2.1.1 Verify Communication

The Touch Prime Ultrasound System is able to verify DICOM communication with a remote system, i.e., PACS and Storage Commitment, RIS, MPPS or a remote printer. This is done at the operator's request. This function is used for network diagnostic purposes.

2.1.2 Store Image and Reports

The Touch Prime Ultrasound System is able to store ultrasound images, structure reports and clips on a remote system. This is done at the operator's request. The remote destination (PACS) can be set in the DICOM setup. Image data to be transferred are instances of the Ultrasound Image Storage or Ultrasound Multi-frame Image Storage SOP class. The structure reports of Vascular and Generic data to be transferred are instances of the Comprehensive SR SOP class.

Graphics, both within the ultrasound image and surrounding the image, are transferred as burned-in graphics, so it is the operator's responsibility to put in or leave out the desired graphics before transferring the image.

2.1.3 Storage Commitment

The Touch Prime Ultrasound System is capable of supporting a remote Storage Commitment System. The user can enable or disable Storage Commitment. If Storage Commitment is configured, the user can choose to configure a separate SCP for Storage Commitment or can configure the Storage SCP to act also as Storage Commitment SCP. The Storage Commitment SCP can be set as part of the PACS Configuration in the DICOM setup.

If Storage Commitment SCP is configured, the Storage SCU issues a Commit request for the images and reports transferred to the Storage SCP by the Storage AE, and the Commitment SCP issues an acknowledgement to the Storage AE on successful commitment of the images and reports.

2.1.4 Retrieve HIS/RIS Worklist

The Touch Prime Ultrasound System is able to retrieve the ultrasound modality Worklist from a RIS. This is done both at the operator's request and automatically at a specified time interval. From the received list, a selection of one Worklist item can be made, i.e. the examination to be performed. The data received from the RIS consists of patient demographic data and procedure step information.

2.1.5 Performed Procedure Step

The Touch Prime Ultrasound System is able to notify a remote Modality Performed Procedure Step System (MPPS) on the procedures performed. The remote system can be configured in the DICOM Setup. It initiates an association with the MPPS to notify the Start, Stop or Discontinued Examinations / Procedures.

2.1.6 Print Image

The Touch Prime Ultrasound System is able to print ultrasound images on a remote grayscale or color printer. The images are sent to the printer at the operator's request or at the end of the current examination.

2.1.7 Query/Retrieve

The Touch Prime Ultrasound System is able to query and retrieve Patient Root Information Model. It issues a set of key attributes to query at the level of patient, study, series or composite object instance. It supplies itself as a move destination for retrieval, and handling the incoming C-STORE sub-operations. It may cancel the ongoing C-FIND or C-MOVE service at any time per user's request.

2.2 Functional Definition of Application Entity

The AE acts as a Service Class User (SCU) for the following SOP classes:

- Verification.
- Storage (PACS and Storage Commitment).
- Basic worklist management.
- Modality Performed Procedure Step (MPPS).
- Basic grayscale print management.
- Basic color print management.
- Patient root query/retrieve information model – FIND

- Patient root query/retrieve information model - MOVE
- Study root query/retrieve information model - MOVE

2.2.1 Real World Activity: Verification

The Verification AE supports as an SCU the following functions:

- Negotiates and establishes DICOM association with remote AEs, i.e., PACS, Storage Commitment Server, RIS, MPPS, and remote printer.
- Verifies communication to a remote AE by issuing an echo request.
- Releases the association with a remote AE.
- Notifies the operator of the communication status.

2.2.2 Real World Activity: Storage

The Storage AE supports as an SCU the following functions:

- Negotiates and establishes DICOM association with a remote PACS (Storage SCP).
- Sends DICOM Information Objects (US, US-mf & Reports) to the remote SCP.
- Requests Commit (Storage Commitment SCP) for the images and reports transferred to the remote PACS.
- Awaits the acknowledgement from the Storage Commitment SCP on successful “Commit” of the images and reports
- Updates the status of the transferred images and reports as “Committed” in the Local Archive (Modality).
- Provides the option of enabling or disabling Storage Commitment.
- The user can choose the PACS SCP to act as both Storage Commitment and PACS or configure a separate Server to act as a Storage Commitment server.
- Releases the association with the remote SCP.
- Notifies the operator of the communication status.

2.2.3 Real World Activity: Request Worklist

The Request Worklist AE supports as an SCU the following functions:

- Negotiates and establishes DICOM association with a remote RIS system (SCP).

- Queries for patient and procedure step information using the Modality Worklist Information Model.
- Releases the association with a remote RIS system.
- Shows the received worklist information to the operator.

2.2.4 Real World Activity: Performed Procedure Step

The Modality Performed Procedure Step AE supports as an SCU the following functions:

- Negotiates and establishes DICOM association with a remote MPPS system (SCP).
- Notifies the MPPS on the Start, Stop and Discontinue of an examination.
- Updates the status of the examination in the local patient database.
- Releases the association with the remote MPPS system.

2.2.5 Real World Activity: Print

The Print AE supports as an SCU the following functions:

- Negotiates and establishes DICOM association with a remote printer (SCP).
- Creates a Film Session.
- Creates one Film Box.
- Sets (updates) one or more Grayscale or Color Image Box.
- Prints (action) a Film Box, i.e., prints one copy of a single film of the film session.
- Releases the association with a remote printer.

2.2.6 Real World Activity: Query/Retrieve

The Q/R AE supports as an SCU the following baseline functions:

- Negotiates and establishes DICOM association with a remote PACS (Storage SCP).
- Issues C-FIND requests that the SCP perform a match of all keys specified in the Identifier of the request down to a Q/R level specified in the request.
- Interpretes various responses for C-FIND command including Pending, Success, Failed, or Refused.
- May cancel the C-FIND service at any time during the process.

- Issues C-MOVE requests providing unique key attribute for each level above the Q/R level based on the query results, with move destination set as itself.
- Handles incoming C-STORE requests as SCP role of C-STORE sub-operations.
- Interpretes various responses for C-MOVE command including Success, Failure or Refused.
- May cancel the C-MOVE service at any time during the process.

2.3 Sequencing of Real-World Activities

Not applicable.

3 AE Specifications

3.1 AE Specification

The Application Entity provides Standard Conformance to the following DICOM SOP classes as an SCU:

Table 3.1-1: Supported SOP Classes by the AE as SCU

SOP Class Name	SOP Class UID
Verification	1.2.840.10008.1.1
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
Storage Commitment Push Model	1.2.840.10008.1.20.1
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31
Basic Grayscale Print Management Meta	1.2.840.10008.5.1.1.9
Basic Color Print Management Meta	1.2.840.10008.5.1.1.18
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3
Comprehensive Structured Report Storage	1.2.840.10008.5.1.4.1.1.88.33
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2

The Application Entity does not support DICOM SOP classes as an SCP, except for storage as sub-operation of Q/R retrieval.

3.1.1 Association Establishment Policies

3.1.1.1 General

The Application Context Name that is always proposed is:

DICOM 3.0 Application Context: "1.2.840.10008.3.1.1.1"

The maximum Protocol Data Unit (PDU) size for PDUs sent to the Touch Prime Ultrasound System is fixed at 32Kb (64Kb for verification). The minimum PDU size accepted for sending

from the scanner is 512 bytes. In case the receiving SCP claims PDUs smaller than 512 bytes, then 512 bytes is used.

3.1.1.2 Number of Associations

The AE will attempt only one association establishment at a time.

3.1.1.3 Asynchronous Nature

The AE does not support asynchronous mode.

3.1.1.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Implementation Class UID: 1.2.208.154.1

Implementation Version Name: BKM DICOM 3.3

3.1.2 Association Initiation Policy

The AE initiates associations for the following activities:

- Verify Communication. See section 3.1.2.1.
- Store Image/Report. See section 3.1.2.2.
- Storage Commitment. See section 3.1.2.3
- Request Worklist. See section 3.1.2.4
- Performed Procedure Step. See section 3.1.2.5
- Print Image. See section 3.1.2.6
- Patient Root Query. See section 3.1.2.7
- Q/R Retrieve. See section 3.1.2.8

3.1.2.1 Real-World Activity: "Echo Operation"

3.1.2.1.1 Associated Real-World Activity

The Associated Real-World Activity is the attempt to verify communications with a remote AE. This occurs when the operator selects the Echo function from the dialog box of the

DICOM setup. In the event that the remote AE does not respond for some reason, the operations will time out after 120 seconds and the association will be released.

3.1.2.1.2 Proposed Presentation Contexts

Table 3.1-2: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1
DICOM Explicit VR Big Endian Transfer Syntax	1.2.840.10008.1.2.2

Table 3.1-3: Proposed Presentation Contexts for Sending Verification Request

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Verification	1.2.840.10008.1.1	See Table 3.1-2	SCU / SCP	None

3.1.2.1.2.1 SOP-Specific Conformance to Verification SOP Class

The AE provides standard conformance. Extended negotiation is not supported.

Verification Service Class is a feature used for network diagnostic purposes. Association is released upon receipt of each C-ECHO confirmation.

3.1.2.2 Real World Activity: Archive Image(s) / Report

3.1.2.2.1 Associated Real-World Activity

The Touch Prime Ultrasound System operator sends a request for storage of an image or a Multi-frame image or a structured report to a remote system. The image or Multi-frame image or structured report is transferred to the remote system.

The remote system is one of the DICOM system settings. After the transfer, the association is released. In the event that the remote system does not respond for some reason, the operations will time out after 120 seconds and the association will be released.

In the event of failure to transfer the Image/Report to the Remote device, the modality tries to re-send the images/reports every 30 seconds, till the configured Max retry attempts are reached. Once the maximum retry limit is reached, the modality stops retrying and notifies the operator on the failed transfer.

The images/reports can either be sent to the Remote system one by one, or as a batch of images/reports together.

3.1.2.2.2 Proposed Presentation Contexts

Table 3.1-4: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1
DICOM Implicit VR Big Endian Transfer Syntax	1.2.830.10008.1.2.2
JPEG Baseline (Process 1): Default Transfer Syntax for Lossy JPEG 8 Bit Image Compression	1.2.840.10008.1.2.4.50
JPEG Lossless (Processes 14): Default Transfer Syntax for Lossless JPEG Image Compression	1.2.840.10008.1.2.4.70

Table 3.1-5: Proposed Presentation Contexts for Image Storage

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	1.2.840.10008.1.2 (default) or 1.2.840.10008.1.2.1 or 1.2.830.10008.1.2.2 See Table 3.1-4	SCU	None
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	1.2.840.10008.1.2.4.50 (default) or 1.2.840.10008.1.2 or 1.2.840.10008.1.2.1 or 1.2.830.10008.1.2.2 See Table 3.1-4	SCU	None

Comprehensive Structured Report Storage	1.2.840.10008.5.1.4.1.1.88.33	1.2.840.10008.1.2 (default) or 1.2.840.10008.1.2.1 or 1.2.830.10008.1.2.2 See Table 3.1-4	SCU	None
---	-------------------------------	---	-----	------

Note: Transfer Syntax for Multi-frame Image Storage or Image Storage can only be changed by trained service personnel. Default for Multi-frame Image Storage is JPEG Baseline (1.2.840.10008.1.2.4.50), and for Image Storage is Implicit VR Little Endian (1.2.840.10008.1.2)

3.1.2.2.2.1 SOP-Specific Conformance to Storage SOP Class

The AE provides standard conformance. Extended negotiation is not supported.

A detailed overview of the applied US Image IOD is given in appendix 9.1.

If a RIS connection is present, Patient and Study related information is retrieved by the AE from the RIS via the Worklist and written in the image headers of the images and Multi-frames to be stored.

The UIDs – (Study Instance UID, Series Instance UID and SOP Instance UID) in the images are generated when the related Study, Series and Image are created. This means that two storages/transfers of the same image will have the same UIDs. The Study Instance UID will be retrieved from the RIS if it is present in the Worklist.

In the following cases, the images and Multi-frames will be resent until the transmission succeeded or the user cancels the jobs:

- If the AE is unable to open an association with the selected destination AE.
- If the Abstract Syntax for an image is not supported by the receiving AE.
- If a failed or refused response to the C-STORE operation is received.

The following are the status codes that are more specifically processed when receiving messages from the Storage SCP equipment:

Table 3.1-6: Storage Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success		0000	
	<i>All others</i>		Association is terminated; Transfer will be retried until aborted by user.

3.1.2.3 Storage Commitment

3.1.2.3.1 Associated Real-World Activity

The Storage AE requests the Storage Commitment of the Storage SOP Classes if a remote AE is configured as Storage Commitment Server (SCP) through the PACS listed in Table 3.1-8 (Storage SCP).

The remote system is one of the DICOM system settings. The Storage commitment acknowledgement can be received in the association requesting the commitment. The Storage AE waits on a configurable listener port, for incoming associations for the Storage Commitment SCP reporting a successful Storage Commitment.

3.1.2.3.2 Proposed Presentation Contexts

The AE will include the presentation context for the Storage Commitment Push Model.

Table 3.1-7: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 3.1-8: Proposed Presentation Contexts for Storage Commitment

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Storage Commitment Push Model	1.2.840.10008.1.20.1	See Table 3.1-7	SCU	None

3.1.2.3.2.1 SOP-Specific Conformance to Storage Commitment SOP Class

The AE provides standard conformance. Extended negotiation is not supported.

A detailed overview of the applied Storage Commitment IOD is given in appendix 9.6.

The Storage AE will send an N-ACTION request to the Storage commitment through the PACS and close the existing association. Then it will wait for a successful event report from the SCP, after successful transfer of an Image or a batch of Images to the PACS. The SCP shall open a new association with the Storage AE and send an N-EVENT-REPORT (Storage Commitment Response) and then the association shall be released by the SCP.

If the report is not received within the applicable time limit for the transaction UID, that specific Storage Commitment will be considered as a failure and the transaction UID is considered as invalid. Whenever the Storage AE tries to resend documents to the PACS, it also resends the Storage Commit requests.

The Storage AE doesn't send the optional image attributes with the Storage Commitment request. The Storage Commitment Status for the various jobs will be stored in the Patient database. The successful / pending / failed Storage commitment status will be updated in the DICOM Status window.

The following table shows the status of the N-ACTION request.

Table 3.1-9: Storage Commitment Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success		0000	
	<i>All others</i>		Association is terminated. The status will be updated as pending till the time out. Once the time out occurs the status will be updated as failure.

The N-EVENT-REPORT transports the status of the Storage Commitment.

3.1.2.4 Retrieve HIS/RIS Worklist

3.1.2.4.1 Associated Real-World Activity

This function can be triggered at the operator’s request or automatically when the Worklist Window is opened. An association is set up to the pre-configured remote system, the RIS. After receiving the Worklist, the association is released. In the event that the remote system does not respond for some reason, the operations will time out after 90 seconds and the association will be released.

3.1.2.4.2 Proposed Presentation Contexts

The AE will include the following presentation contexts:

Table 3.1-10: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1

Table 3.1-11: Proposed Presentation Contexts for Request for Modality Worklist

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	See Table 3.1-10	SCU	None

3.1.2.4.2.1 SOP Specific Conformance to Modality Worklist Management

The AE provides standard conformance. Extended negotiation is not supported.

An overview of all requested Matching and Return keys with additional attribute information is given in appendix 9.2. The matching type (Single Value, Wild Card Matching or Range Matching) is also specified.

The user can chose between four date/time queries and one patient query.

The four date/time queries are: Today only, +/- 12 hours, +/- 24 Hours and +/- 3 Days.

The patient query is made from Patient Name (only Last is used), Patient ID, Accession Number and Requested Procedure ID. The user can enter data in one or more of these fields to query for patients i.e. the user can enter data in the Accession Number field only to query on Accession Number only.

The system will expect the extended character set in the worklist (used on the RIS) to match the extended character set on the scanner. The user will be notified if the two extended character sets do not match. The extended character set used on the scanner depends on the selected language. See section 7 for a list of language and extended character sets.

The following are the status codes that are more specifically processed when messages are received from the Modality Worklist SCP equipment:

Table 3.1-12: Modality Worklist Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Refused	Out of Resources	A700	Process terminated. No patient data received.
Failed	Identifier does not match SOP Class	A900	Process terminated. No patient data received.
	Unable to process	Cxxx	Process terminated. No patient data received.
Cancel	Matching terminated due to Cancel request	FE00	Process terminated. No patient data received.
Success	Matching is complete – No final identifier is supplied	0000	
Pending	Matches are continuing – Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys.	FF00	The process of receiving matches continues.

	Matches are continuing – Warnings that one or more Option Keys were not supported for existence and/or matching for this identifier.	FF01	The process of receiving matches continues without any warnings or errors.
	Unknown	None of above	Process terminated. No patient data received.

If the response of a query to the RIS/Worklist-server somehow fails, the user will see an empty list of patients and examinations. The operator can retry by pushing an update key or can enter the patient and examination information manually.

3.1.2.5 Performed Procedure(s)

3.1.2.5.1 Associated Real-World Activity

An association is created with the specific remote Modality Performed Procedure Step (MPPS) System when the user selects the patient and on OK in the Patient Entry window or triggers a new examination for the existing patient. The related MPPS SOP instance is created at once.

The status message is sent to the MPPS during Start, Stop or Discontinued Examination. After every message transfer the association is released.

3.1.2.5.2 Proposed Presentation Contexts

The AE will include the following presentation contexts for notification of the performed procedures:

Table 3.1-13: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 3.1-14: Proposed Presentation Contexts for MPPS

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	See Table 3.1-13	SCU	None

3.1.2.5.2.1 SOP Specific Conformance to Modality Performed Procedure Step SOP Classes

The AE provides standard conformance. Extended negotiation is not supported.

A detailed overview of the applied Modality Performed Procedure Step IOD is given in appendix 9.5.

The N-CREATE Service Element is used to create the MPPS instance. The N-SET Service Element is used to indicate the end of the MPPS.

Table 3.1-15: MPPS Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success		0000	
	<i>All others</i>		Association is terminated; the status is updated as failure.

3.1.2.6 Print Image(s)

3.1.2.6.1 Associated Real-World Activity

An association is initiated with the named DICOM printer when the operator requests the image to be printed. After the printing is finished, the association is released. In the event that the printer does not respond for some reason, the operations will time out after 120 seconds and the association will be released.

3.1.2.6.2 Proposed Presentation Contexts

The AE will include the following presentation contexts for printing an image as a grayscale image:

Table 3.1-16: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2

Table 3.1-17: Proposed Presentation Contexts for Grayscale Print Management

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Basic Grayscale Print Management Meta SOP Class	1.2.840.10008.5.1.1.9	See Table 3.1-16	SCU	None

Table 3.1-18: Proposed Presentation Contexts for Color Print Management

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Basic Color Print Management Meta SOP Class	1.2.840.10008.5.1.1.18	See Table 3.1-16	SCU	None

3.1.2.6.2.1 SOP Specific Conformance to Print SOP Classes

The AE provides standard conformance. Extended negotiation is not supported.

The N-CREATE Service Element is used for the Basic Film Session and Basic Film Box. The N-SET Service Element is used for Basic Grayscale Image Box and Basic Color Image Box. The N-ACTION Service Element is used for Basic Film Box to print the image. The N-EVENT-REPORT Service Element is used to report the changes of printer status in an asynchronous way.

The color/monochrome configuration must be correctly set for the printer. A color printer that is set up as a monochrome printer (or vice versa) will not produce any output. A printer that supports both color and monochrome must be installed as a color printer.

The implementation does make use of certain User Optional attributes that are mandatory for the DICOM printer. See appendix 9.3 and 9.4 for details.

The following are the status codes that are more specifically processed when messages are received from the Print SCP equipment:

Table 3.1-19: Create Basic Film Session Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success	Film session successfully created	0000	
	<i>All others</i>		Association is terminated; Creation will be retried until aborted by user.

Table 3.1-20: Create Basic Film Box Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success	Film box successfully created	0000	
	<i>All others</i>		Association is terminated; Creation will be retried until aborted by user.

Table 3.1-21: Set Basic Grayscale and Color Image Box Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success	Image successfully stored in image box	0000	
	<i>All others</i>		Association is terminated; Update (set) will be retried until aborted by user.

Table 3.1-22: Print Basic Film Box Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Success	Film accepted for printing	0000	
	<i>All others</i>		Association is terminated; Print job will be retried until aborted by user.

Table 3.1-23: Print SOP Class Status Codes

Service Status	Further Meaning	Event type ID	Application behavior when receiving Status Codes
NORMAL	Given jobs are being printed	1	
	<i>All others</i>		Association is terminated; Print job will be retried until aborted by user.

3.1.2.7 Patient Root Query

3.1.2.7.1 Associated Real-World Activity

This function is triggered at the operator’s request after entering query values for a given set of key attributes at different level of patient, study and series. All matching patients are received and displayed. Selecting a patient triggers query of studies/series for the selected patient.

An association is set up to the pre-configured default DICOM storage device. After query is complete, the association is released.

3.1.2.7.2 Proposed Presentation Contexts

The AE will include the following presentation contexts:

Table 3.1-24: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1
DICOM Explicit VR Big Endian Transfer Syntax	1.2.840.10008.1.2.2

Table 3.1-25: Proposed Presentation Contexts for Request for Patient Root Q/R Information Model - FIND

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1	See Table 3.1-1024	SCU	None

3.1.2.7.2.1 SOP Specific Conformance to Patient Root Q/R Information Model - FIND

The AE provides baseline behavior of the C-FIND SCU. Extended negotiation is not supported.

An overview of all presented Key Attributes with additional attribute information is given in appendix 9.8. The matching type (Single Value, List of UID, Universal Matching, Wild Card Matching, Range Matching or Sequence Matching) is also specified. Particularly, patient name and patient ID are wild card matching.

The user can enter data for one or more of the query key attributes. Note some the key attributes may not be visible for some vendors even if they are supported in the software.

The system will expect the extended character set in the key attributes to match the extended character set on the scanner. See section 7 for a list of language and extended character sets.

The following are the specific status codes which might be returned in a C-FIND response.

Table 3.1-26: Patient Root Q/R Information Model - FIND Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Failed	Refused: Out of Resources	A700	Process terminated. No data received.
	Identifier does not match SOP Class	A900	Process terminated. No data received.

	Unable to process	Cxxx	Process terminated. No data received.
Cancel	Matching terminated due to Cancel request	FE00	Process terminated. Stopped receiving remaining data.
Success	Matching is complete – No final identifier is supplied	0000	Complete process.
Pending	Matches are continuing – Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys.	FF00	The process of receiving matches continues.
	Matches are continuing – Warnings that one or more Option Keys were not supported for existence and/or matching for this identifier.	FF01	The process of receiving matches continues without any warnings or errors.

If the response of a query to the Patient Root Q/R Information Model somehow fails or no matching records returned, the user will see an empty list of patients and examinations. The operator can modify the query key values then retry the search.

3.1.2.8 Q/R Retrieve

3.1.2.8.1 Associated Real-World Activity

This function is triggered when a queried study is selected. All the ultrasound series and instances belonging to the selected study will be retrieved.

Both Study Root MOVE and Patient Root Move are supported. It might be configured to try Study Root MOVE first before Patient Root MOVE. Otherwise, Patient Root MOVE is always executed.

An association is set up to the pre-configured default DICOM storage device. After C-MOVE command is complete, the association is released. The subsequent C-STORE sub-operation shall occur on a different association where the local AE serves as an SCP of the storage service class.

3.1.2.8.2 Proposed Presentation Contexts

The AE will include the following presentation contexts:

Table 3.1-27: Transfer Syntax

Name List	UID List
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1
DICOM Explicit VR Big Endian Transfer Syntax	1.2.840.10008.1.2.2

Table 3.1-28: Proposed Presentation Contexts for Request for Q/R Information Model - MOVE

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2	See Table 3.1-1027	SCU	None
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2	See Table 3.1-1027	SCU	None

3.1.2.8.2.1 SOP Specific Conformance to Patient Root or Study Root Q/R Information Model - MOVE

The AE provides baseline behavior of the C-MOVE SCU for either Patient Root or Study Root Q/R Information Model. Extended negotiation is not supported.

The Patient Root C-MOVE is the default method. It can be configured to first try Study Root C-MOVE, and if it is not supported by Q/R SCP, it falls through to Patient Root C-MOVE. When Study Root C-MOVE is issued, the AE is listening on the local SCU port, default 104, for incoming association of C-STORE sub-operations. If Patient Root C-MOVE is used, the AE is listening on a configurable local port, default 7878, for incoming C-STORE association. If the port is the same as local SCU port, the destination AE title of C-MOVE is set to be same as local AE title. Otherwise, the destination AE title is the local AE title suffixed with “_SCP”, e.g. if the local AE title is “ABC” then the target AE title would be “ABC_SCP”.

A detailed overview of the applied C-MOVE IOD is given in appendix 9.9.

The following are the specific status codes which might be returned in a C-MOVE response.

Table 3.1-29: Q/R Information Model – MOVE Status Codes

Service Status	Further Meaning	Status Codes	Application behavior when receiving Status Codes
Failure	Refused: Out of Resources – Unable to calculate number of matches	A701	Process terminated. No data received.
	Refused: Out of Resources – Unable to perform sub-operations	A702	Process terminated. No data received.
	Refused: Move Destination unknown	A801	Process terminated. No data received.
	Identifier does not match SOP Class	A900	Process terminated. No data received.
	Unable to process	Cxxx	Process terminated. No data received.
Cancel	Sub-operation terminated due to Cancel Indication	FE00	Process terminated. Stopped receiving remaining data.
Warning	Sub-operation Complete – One or more Failures	B000	Complete process.
Success	Sub-operation Complete – No Failures	0000	Complete process.
Pending	Sub-operations are continuing	FF00	Continue process.

If the response to the C-MOVE somehow fails or no matching composite object instances returned, the user will see an empty list of thumbnails.

3.1.3 Association Acceptance Policy

The AE does not accept associations proposed by another application entity.

4 Communication Profiles

4.1 TCP/IP Stack

The AE provides DICOM 3.0 TCP/IP Network Communication Support as defined in DICOM PS 3.8.

The TCP/IP stack is inherited from the underlying Microsoft Windows®-based operating system.

4.1.1 Physical Media Support

The system is indifferent to the physical medium; it inherits the medium from the Windows® Operating System.

5 Extensions / Specializations / Privatizations

Not applicable.

6 Configuration

The system is configured in the DICOM setup menu. Only an operator with the correct level of knowledge should change the configuration. The Configuration menu is intended to be used during installation, by a service engineer only.

6.1 AE Title/Presentation Address Mapping

The Local AE Title is configurable. A service engineer must configure it during installation.

6.2 Configurable Parameters

6.2.1 Local AE

The following fields are configurable:

- Local AE Title.

The following fields can be configured through Windows Network Setup dialog:

- Local IP Address.
- Local IP Net mask.

6.2.2 Remote AE

The following fields are configurable for every remote DICOM AE:

- Remote AE Title.
- Remote IP Address or hostname.
- Responding TCP/IP Port.

6.2.3 Storage

The default transfer syntax for Image Storage is Implicit VR Little Endian. The transfer syntax can be changed to Explicit VR Little Endian by Carestream trained service personnel.

The default transfer syntax for Multi-frame images (clips) is JPEG. The transfer syntax can be changed to Implicit Little Endian by Carestream trained service personnel.

The default transfer syntax for Reports Storage is Explicit VR Little Endian. The transfer syntax can be changed to Implicit VR Little Endian by Carestream trained service personnel.

6.2.4 Storage Commitment

The Storage Commitment support for the Modality can be configured. By default Storage Commitment support is disabled. The Storage Commitment support can be enabled or disabled & either the Storage SCP can be configured as Storage Commitment SCP or a new Storage Commitment SCP can be configured. For the Storage commitment SCU, a listener port can be configured. Only Carestream trained service personnel can do these changes.

6.2.5 Printing

The color/monochrome configuration must be correctly set for the printer. A color printer that is set up as a monochrome printer (or vice versa) will not produce any output. A printer that supports both color and monochrome must be installed as a color printer.

6.2.6 Query/Retrieve

The port number for the C-MOVE destination AE to listen on for incoming C-STORE sub-operations is configurable by service engineer. The default port number is 7878. Note that if this port is equal to the local AE port, the AE title for the destination AE is set the same as local AE. Otherwise, its title is the local AE title suffixed with “_SCP”.

Patient Root Q/R Information Model is the default to use for C-MOVE. It can be configured to try Study Root first by service engineer. But be aware in this case, the listening port specified above is obsolete, as the local AE port now takes precedence.

7 Support of Extended Character Sets

The system supports the following character sets depending upon the language selected on the scanner:

Language	DICOM Character Set
BULGARIAN	ISO_IR 144
CHINESE	ISO_IR 192
CZECH	ISO_IR 101
<i>DEFAULT¹</i>	ISO_IR 100 / ISO_IR 6 ²
ESTONIAN	ISO_IR 110
GREEK	ISO_IR 126
HUNGARIAN	ISO_IR 101
LATVIAN	ISO_IR 110
LITHUANIAN	ISO_IR 110
POLISH	ISO_IR 101
ROMANIAN	ISO_IR 101
RUSSIAN	ISO_IR 144
SLOVAK	ISO_IR 101

8 Security Profiles

The product does not conform to any defined DICOM Security Profiles.

¹ Default refers to the following languages:

Danish, Dutch, English, Finnish, French, German, Icelandic, Italian, Norwegian, Portuguese, Spanish, Swedish

² ISO_IR 6 if no extended characters are used.

It is assumed that the product is used within a secured environment.

It is assumed that a secured environment includes at a minimum:

1. Firewall or router protections to ensure that only approved external hosts have network access to the product.
2. Firewall or router protections to ensure that the product only has network access to approved external hosts and services.
3. Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN)).

9 Appendix: Overview of the Applied Information Object Definitions

Definitions

Usage Specification M = Mandatory C = Conditional U = User Option

Matching Key U = Unique R = Required O = Optional

Type/Return Key 1 = Mandatory 2 = Mandatory, may be empty 3 = Optional

 1C = Conditional 2C = Conditional

Device usage A = Always AE = Always maybe empty N = Not used

 No entry = standard behaviour

9.1 Overview of the Applied Ultrasound (US) Storage SOP Classes

The modules selected from the IOD module table of DICOM 3.0 are given in the table below.

Table 9.1-1: Applied Modules in the US Image IOD

Information Entity	Module	Usage	Details
Patient	Patient	M	See Table 9.1-3
Study	General Study	M	See Table 9.1-4
	Patient Study	U	See Table 9.1-5
Series	General Series	M	See Table 9.1-6
Equipment	Equipment	M	See Table 9.1-7
Image	General Image	M	See Table 9.1-8
	Image Pixel	M	See Table 9.1-9
	Contrast/Bolus	C	See Table 9.1-10
	US Region Calibration	U	See Table 9.1-14
	US Image	M	See Table 9.1-15
	SOP Common	M	See Table 9.1-16

Table 9.1-2: Applied Modules in the US Multi-frame Image IOD

Information Entity	Module	Usage	Details
Patient	Patient	M	See Table 9.1-3
Study	General Study	M	See Table 9.1-4
	Patient Study	U	See Table 9.1-5
Series	General Series	M	See Table 9.1-6
Equipment	Equipment	M	See Table 9.1-7
Image	General Image	M	See Table 9.1-8
	Image Pixel	M	See Table 9.1-9
	Contrast/Bolus	C	See Table 9.1-10
	Cine	M	See Table 9.1-11
	Multi-frame	M	See Table 9.1-12
	Frame Pointers	U	See Table 9.1-13
	US Image	M	See Table 9.1-15
SOP Common	M	See Table 9.1-16	

The details of these applied modules are given in the tables below. The list of possible values is given, if applicable. If an attribute may be present conditionally/optionally or may contain a zero length value, this is also indicated.

Table 9.1-3: Patient Module

Attribute Name	Tag	Type	Note
Patient's Name	(0010,0010)	2A	Received from RIS or entered by user
Patient ID	(0010,0020)	2A	Received from RIS or entered by user
Patient's Birth Date	(0010,0030)	2	Received from RIS or entered by user
Patient's Sex	(0010,0040)	2	Received from RIS or entered by user
Other Patient ID's	(0010,1000)	2	Received from RIS or entered by user

Table 9.1-4: General Study Module

Attribute Name	Tag	Type	Note
Study Instance UID	(0020,000D)	1	Generated at creation of the Study or received from RIS
Study Date	(0008,0020)	2A	Generated at creation of the Study
Study Time	(0008,0030)	1	Generated at creation of the Study
Referring Physician's Name	(0008,0090)	2	Received from RIS or entered by user

Study ID	(0020,0010)	2	Auto-generated from Study Date and Study Time if not entered by user
Accession Number	(0008,0050)	2	Received from RIS or entered by user
Study Description	(0008,1030)	3AE	Received from RIS (copied from Scheduled Procedure Step Description (0040,0007)) or entered by the user
Name of Physician(s) Reading Study	(0008,1060)	3AE	Received from RIS (copied from Names of Intended Recipients of Results (0040,1010)) or entered by the user

Table 9.1-5: Patient Study Module

Attribute Name	Tag	Type	Note
Admitting Diagnosis Description	(0008,1080)	3AE	Present if received from RIS or entered by user or empty
Patient's Size	(0010,1020)	3	Length or size of the Patient, in meters. Not present if neither received from RIS nor entered by user
Patient's Weight	(0010,1030)	3	Weight of the Patient, in kilograms. Not present if neither received from RIS nor entered by user
Admission ID	(0038,0010)	3	Present if received from RIS or entered by user or empty

Table 9.1-6: General Series Module

Attribute Name	Tag	Type	Note
Modality	(0008,0060)	1	US or SR
Series Instance UID	(0020,000E)	1	Generated at creation of the Series
Series Number	(0020,0011)	2	Auto generated
Performing Physician's Name	(0008,1050)	3	Entered by user
Protocol Name	(0018,1030)	3	User defined
Operator's Name	(0008,1070)	3	Entered by user
Body Part Examined	(0018,0015)	3	Entered by user
Requested Attribute Sequence	(0040,0275)	3AE	Active only if MPPS & MWL are available.
Requested Procedure ID	(0040,1001)	1C	From MWL
Scheduled Procedure Step ID	(0040,0009)	1C	From MWL

Scheduled Procedure Step Description	(0040,0007)	3	From MWL
Performed Procedure Step ID	(0040,0253)	3	Auto generated – only if MPPS is active
Performed Procedure Step Start Date	(0040,0244)	3	Auto generated– only if MPPS is active
Performed Procedure Step Start Time	(0040,0245)	3	Auto generated– only if MPPS is active
Performed Procedure Step Description	(0040,0254)	3	Entered by user – active only if MPPS is available.
Performed Protocol Code Sequence	(0040,0260)	3	Auto generated – only if MPPS is active

Table 9.1-7: General Equipment Module

Attribute Name	Tag	Type	Note
Manufacturer	(0008,0070)	2A	“Carestream Health”
Institution Name	(0008,0080)	3AE	Institution where the equipment that produced the composite instance is located. Set up by the user.
Institution Address	(0008,0081)	3	Mailing address of the institution where the equipment that produced the composite instances is located. Not used
Station Name	(0008,1010)	3	User defined name identifying the machine that produced the composite instances. Not used
Institutional Department Name	(0008,1040)	3	Department in the institution where the equipment that produced the composite instances is located. Not used.
Manufacturer’s Model Name	(0008,1090)	3A	“Touch Ultra Sound System”
Device Serial Number	(0018,1000)	3A	The serial number of the scanner
Software Version(s)	(0018,1020)	3A	Software version of the current software

Table 9.1-8: General Image Module

Attribute Name	Tag	Type	Note
Instance Number	(0020,0013)	2	Images, clips and reports are numbered in chronological order, starting from 1 at the beginning of a new examination
Patient Orientation	(0020,0020)	2C	
Content Date	(0008,0023)	2C	Date when the image was captured.

			Note: This Attribute was formerly known as Image Date
Content Time	(0008,0033)	2C	Time when the image was captured. Note: This Attribute was formerly known as Image Time
Image Type	(0008,0008)	3	
Burned In Annotation	(0028,0301)	3A	"YES"
Lossy Image Compression	(0028,2110)	3	01H for Multi-frame Image with JPEG Transfer Syntax. Otherwise: Not present.

Table 9.1-9: Image Pixel Module

Attribute Name	Tag	Type	Note
Samples per Pixel	(0028,0002)	1	0003H
Photometric Interpretation	(0028,0004)	1	Multi-frame with JPEG Transfer Syntax: "YBR_FULL_422" Otherwise: "RGB"
Rows	(0028,0010)	1	Image height
Columns	(0028,0011)	1	Image width
Bits Allocated	(0028,0100)	1	0008H
Bits Stored	(0028,0101)	1	0008H
High Bit	(0028,0102)	1	0007H
Pixel Representation	(0028,0103)	1	0000H (= unsigned integer)
Planar Configuration	(0028,0006)	1C	1 (color-by-plane)
Pixel Data	(7FE0,0010)	1	

Table 9.1-10: Contrast/Bolus Module

Attribute Name	Tag	Type	Note
Contrast/Bolus Agent	(0018,0010)	2	Present and empty if Contrast Harmonic is activated. Otherwise not present.

Table 9.1-11: Cine Module

Attribute Name	Tag	Type	Note
Preferred Playback Sequencing	(0018,1244)	3	0 (looping playback)
Frame Time	(0018,1063)	1C	1000 / Frame rate of Multi-frame in Hz (Nominal time in msec per individual

			frame)
Recommended Display Frame Rate	(0008,2144)	3	Frame rate of Multi-frame in Hz (Recommended rate (frames/sec) for display of multi-frame sequence)
Cine Rate	(0018,0040)	3	Frame rate of Multi-frame in Hz (Number of frames per seconds)
Frame Delay	(0018,1066)	3	000H (Time in msec to start of first frame)

Table 9.1-12: Multi-frame Module

Attribute Name	Tag	Type	Note
Number of Frames	(0028,0008)	1	Number of frames in Multi-frame
Frame Increment Pointer	(0028,0009)	1	Set to Frame Time (0018,1063) from Table 9.1-11: Cine Module

Table 9.1-13: Frame Pointers Module

Attribute Name	Tag	Type	Note
Representative Frame Number	(0028,6010)	3	1 (Number of frame selected for use as icon)

Table 9.1-14: US Region Calibration Module

Attribute Name	Tag	Type	Note
Sequence of Ultrasound Regions	(0018,6011)	1	A sequence for each B-mode view. One sequence for each D/M/CW view.
>Region Location Min x0	(0018,6018)	1	
>Region Location Min y0	(0018,601A)	1	
>Region Location Max x1	(0018,601C)	1	
>Region Location Max y1	(0018,601E)	1	
>Physical Units X Direction	(0018,6024)	1	B-mode: 03H (cm) D-mode: 04H (sec) M-mode: 04H (sec) CW-mode: 04H (sec)
>Physical Units Y Direction	(0018,6026)	1	B-mode: 03H (cm) D-mode: 07H (cm/sec) M-mode: 03H (cm) CW-mode: 07H (cm/sec)
>Physical Delta X	(0018,602C)	1	

>Physical Delta Y	(0018,602E)	1	
>Reference Pixel x0	(0018,6020)	3	
>Reference Pixel y0	(0018,6022)	3	
>Ref. Pixel Physical Value X	(0018,6028)	3	0
>Ref. Pixel Physical Value Y	(0018,602A)	3	0
>Region Spatial Format	(0018,6012)	1	B-mode: 01H (2D) D-mode: 03H (Spectral) M-mode: 02H (M-Mode) CW-mode: 03H (Spectral)
>Region Data Type	(0018,6014)	1	B-mode: 01H (Tissue) D-mode: 03H (PW Spectral Doppler) M-mode: 01H (Tissue) CW-mode: 04H (CW Spectral Doppler)
>Region Flags	(0018,6016)	1	0
>Doppler Correction Angle	(0018,6034)	3	B-mode: Not applied. D-mode: Angle Correction (Degrees) M-mode: Not applied. CW-mode: Not applied

Table 9.1-15: US Image Module

Attribute Name	Tag	Type	Note
Samples per Pixel	(0028,0002)	1	See Table 9.1-9 Image Pixel Module
Photometric Interpretation	(0028,0004)	1	See Table 9.1-9 Image Pixel Module
Bits Allocated	(0028,0100)	1	See Table 9.1-9 Image Pixel Module
Bits Stored	(0028,0101)	1	See Table 9.1-9 Image Pixel Module
High Bit	(0028,0102)	1	See Table 9.1-9 Image Pixel Module
Planar Configuration	(0028,0006)	1C	See Table 9.1-9 Image Pixel Module
Pixel Representation	(0028,0103)	1	See Table 9.1-9 Image Pixel Module
Frame Increment Pointer	(0028,0009)	1C	See Table 9.1-9 Image Pixel Module
Image Type	(0008,0008)	2	See Table 9.1-8: General Image Module
Lossy Image Compression	(0028,2110)	1C	See Table 9.1-8: General Image Module
Transducer Data	(0018,5010)	3	Name of transducer e.g. "8801". Only set for single frame images
Mechanical Index	(0018,5022)	3	Only sent for single view images
Bone Thermal Index	(0018,5024)	3	Only sent for single view images. Only sent if TI type is TIB.
Cranial Thermal Index	(0018,5026)	3	Only sent for single view images.

			Only sent if TI type is TIC.
Soft Tissue Thermal Index	(0018,5027)	3	Only sent for single view images. Only sent if TI type is TIS.

Table 9.1-16: SOP Common Module

Attribute Name	Tag	Type	Note
SOP Class UID	(0008,0016)	1	For US Image: "1.2.840.10008.5.1.4.1.1.6.1" For US Multi-Frame Image: "1.2.840.10008.5.1.4.1.1.3.1"
SOP Instance UID	(0008,0018)	1	Generated when image is created
Specific Character Set	(0008,0005)	1C	Set according to selected language on the scanner. See section 7.
Instance Creation Date	(0008,0012)	3	Document creation date
Instance Creation Time	(0008,0013)	3	Document creation time
Time zone Offset From UTC	(0008,0201)	3	
Instance Number	(0020,0013)	3	See Table 9.1-8: General Image Module

9.2 Overview of the Applied Modality Worklist IOD

This section specifies in detail the applied attributes in the C-FIND Service Element of this supported SOP Class.

If an attribute is present conditionally/optionally or if the attribute may contain a zero length value, this is indicated.

The scanner will use predefined DICOM Character Sets depending upon the selected language. See table in section 7.

The user will be warned if the Character Set of the worklist does not match the Character Set of the scanner.

The search filter mentioned is set in the worklist setup for the modality. The queries mentioned are selected in the patient dialog where the worklist is displayed.

Table 9.2-1: Modality Worklist Information Model - FIND SOP Class - C-FIND

Description	Tag	Match	Return	Note
Scheduled Procedure Step				
Scheduled Procedure Step Sequence	(0040,0100)	R	1	Return key
> Scheduled Station AE Title	(0040,0001)	R	1	Single Value Matching is applied; the applied value is the configured AE Title. Only used if search filter is "This System Only"
> Scheduled Procedure Step Start Date	(0040,0002)	R	1	Range Value Matching is applied when using date/time queries.
> Scheduled Procedure Step Start Time	(0040,0003)	R	1	Range Value Matching is applied when using date/time queries.
> Modality	(0008,0060)	R	1	Single Value Matching is applied; the applied value is US. Only used if search filter is "Ultrasound Only"
> Scheduled Procedure Step Description	(0040,0007)	O	1C	Return key
> Scheduled Procedure Step Location	(0040,0011)	O	2	Return key
> Scheduled Protocol Code Sequence	(0040,0008)	O	1C	Return key
> Scheduled Procedure Step ID	(0040,0009)	O	1	Return key
Requested Procedure				
Requested Procedure ID	(0040,1001)	R	1	Return key. Single Value Matching is applied when data has been entered by the user and "Patient" query has been selected.
Requested Procedure Description	(0032,1060)	R	1C	Return Key
Requested Procedure Code Sequence	(0032,1064)	R	1C	Return Key
> Code Value	(0008,0100)	R	1	Return Key
> Coding Scheme Designator	(0008,0102)	R	1	Return Key
> Code Meaning	(0008,0104)	R	2	Return key
Study Instance UID	(0020,000D)	R	1	Return key
Referenced Study	(0008,1110)	R	2	Return Key

Sequence				
> Referenced SOP Class UID	(0008,1150)	R	1	Return Key
> Referenced SOP Instance UID	(0008,1155)	R	1	Return Key
Requested Procedure Priority	(0040,1003)	O	2	Return key
Names of Intended Recipients of Results	(0040,1010)	O	3	Return key
Imaging Service Request				
Accession Number	(0008,0050)	R	1	Return key. Single Value Matching is applied when data has been entered by the user and "Patient" query has been selected.
Referring Physician's Name	(0008,0090)	R	2	Return key
Visit Admission				
Admitting Diagnosis Description	(0008,1080)	O	2	Return key
Admission ID	(0038,0010)	O	3	Return key
Patient Identification				
Patient's Name	(0010,0010)	R	1	Return key. Single Value Matching is applied when data has been entered by the user and "Patient" query has been selected.
Patient ID	(0010,0020)	R	1	Return key. Single Value Matching is applied when data has been entered by the user and "Patient" query has been selected.
Other Patient ID's	(0010,1000)	R	2	Return key.
Patient Demographic				
Patient's Birth Date	(0010,0030)	R	2	Return key
Patient's Sex	(0010,0040)	R	2	Return key
Patient's Weight	(0010,1030)	R	2	Return key
Patient's Size	(0010,1020)	R	2	Return key
Patient's Address	(0010,1040)	O	3	Return key
Current Patient Location	(0038,0300)	O	3	Return key

9.3 Overview of the Applied Basic Grayscale Print Management Meta IOD

9.3.1 Overview of the Applied Basic Film Session IOD

This section specifies in detail the applied attributes in the N-CREATE Service Element of this supported SOP Class.

Table 9.3-1: Basic Film Session Presentation Module

Attribute Name	Tag	Usage	Note
Number of Copies	(2000,0010)	U	1
Print Priority	(2000,0020)	U	HIGH
Medium Type	(2000,0030)	U	PAPER, CLEAR FILM, BLUE FILM
Film Destination	(2000,0040)	U	PROCESSOR
Film Session Label	(2000,0050)	U	"Session <month>-<date>-<year>-<hours>:<minutes>"

9.3.2 Overview of the Applied Basic Film Box IOD

This section specifies in detail the applied attributes in the N-CREATE Service Element of this supported SOP Class.

Table 9.3-2: Basic Film Box Presentation Module

Attribute Name	Tag	Usage	Note
Image Display Format	(2010,0010)	M	STANDARD\C,R (C,R = 1,1 / 1,2 / 2,2 / 2,3 / 3,3 / 3,4 / 3,5 / 4,4 / 4,5 / 4,6 / 5,5 / 5,6 / <Custom>)
Referenced Film Session Sequence	(2010,0500)	M	
> Referenced SOP Class UID	(0008,1150)	M	Appl. value: 1.2.840.10008.5.1.1.1 (Basic Film Session SOP Class)
> Referenced SOP Instance UID	(0008,1155)	M	Appl. value: The SOP Instance UID of the parent film session
Film Orientation	(2010,0040)	U	PORTRAIT / LANDSCAPE
Film Size ID	(2010,0050)	U	See defined terms
Magnification Type	(2010,0060)	U	CUBIC

Film Size ID

The defined terms are:

8INX10IN 10INX12IN

10INX14IN 11INX14IN

14INX14IN 14INX17IN

24CMX24CM 24CMX30CM

and custom size (both CM and IN)

9.3.3 Overview of the Applied Basic Grayscale Image Box IOD

This section specifies in detail the applied attributes in the N-SET Service Element of this supported SOP Class.

Table 9.3-3: Basic Grayscale Image Box Presentation Module

Attribute Name	Tag	Usage	Note
Image Position	(2020,0010)	M	
Polarity	(2020,0020)	U	NORMAL
Basic Grayscale Image Sequence	(2020,0110)	M	
> Samples Per Pixel	(0028,0002)	M	0001H
> Photometric Interpretation	(0028,0004)	M	MONOCHROME2
> Rows	(0028,0010)	M	Image height
> Columns	(0028,0011)	M	Image width
> Pixel Aspect Ratio	(0028,0034)	U	"1\1"
> Bits Allocated	(0028,0100)	M	0008H
> Bits Stored	(0028,0101)	M	0008H
> High Bit	(0028,0102)	M	0007H
> Pixel Representation	(0028,0103)	M	0000H (= unsigned integer)
> Pixel Data	(7FE0,0010)	M	

9.4 Overview of the Applied Basic Color Print Management Meta IOD

9.4.1 Overview of the Applied Basic Film Session IOD

This section specifies in detail the applied attributes in the N-CREATE Service Element of this supported SOP Class.

Table 9.4-1: Basic Film Session Presentation Module

Attribute Name	Tag	Usage	Note
Number of Copies	(2000,0010)	U	1
Print Priority	(2000,0020)	U	HIGH
Medium Type	(2000,0030)	U	PAPER, CLEAR FILM, BLUE FILM
Film Destination	(2000,0040)	U	PROCESSOR
Film Session Label	(2000,0050)	U	“Session <month>-<date>-<year>-<hours>:<minutes>”

9.4.2 Overview of the Applied Basic Film Box IOD

This section specifies in detail the applied attributes in the N-CREATE Service Element of this supported SOP Class.

Table 9.4-2: Basic Film Box Presentation Module

Attribute Name	Tag	Usage	Note
Image Display Format	(2010,0010)	M	STANDARD\C,R (C,R = 1,1 / 1,2 / 2,2 / 2,3 / 3,3 / 3,4 / 3,5 / 4,4 / 4,5 / 4,6 / 5,5 / 5,6 / <Custom>)
Referenced Film Session Sequence	(2010,0500)	M	
> Referenced SOP Class UID	(0008,1150)	M	Appl. value: 1.2.840.10008.5.1.1.1 (Basic Film Session SOP Class)
> Referenced SOP Instance UID	(0008,1155)	M	Appl. value: The SOP Instance UID of the parent film session
Film Orientation	(2010,0040)	U	PORTRAIT / LANDSCAPE
Film Size ID	(2010,0050)	U	See defined terms
Magnification Type	(2010,0060)	U	CUBIC

Film Size ID

The defined terms are:

- 8INX10IN 10INX12IN
- 10INX14IN 11INX14IN
- 14INX14IN 14INX17IN
- 24CMX24CM 24CMX30CM

and custom size (both CM and IN)

9.4.3 Overview of the Applied Basic Color Image Box IOD

This section specifies in detail the applied attributes in the N-SET Service Element of this supported SOP Class.

Table 9.4-3: Basic Color Image Box Presentation Module

Attribute Name	Tag	Usage	Note
Image Position	(2020,0010)	M	
Polarity	(2020,0020)	U	NORMAL
Basic Color Image Sequence	(2020,0111)	M	
> Samples Per Pixel	(0028,0002)	M	3
> Photometric Interpretation	(0028,0004)	M	RGB
> Rows	(0028,0010)	M	Image height
> Columns	(0028,0011)	M	Image width
> Pixel Aspect Ratio	(0028,0034)	U	"1\1"
> Bits Allocated	(0028,0100)	M	0008H
> Bits Stored	(0028,0101)	M	0008H
> High Bit	(0028,0102)	M	0007H
> Pixel Representation	(0028,0103)	M	0000H (= unsigned integer)
> Pixel Data	(7FE0,0010)	M	

9.5 Overview of the Applied Modality Performed Procedure Step IOD

This section specifies in detail the applied attributes in the N-CREATE and N-SET Service Elements of this supported SOP Class.

Table 9.5-1: Performed Procedure Step Relationship Module attributes

Attribute Name	Tag	Req.Type N-CREATE	Req.Type N-SET	Note
Scheduled Step Attributes Sequence	(0040,0270)	1	Not allowed	From MWL
>Study Instance UID	(0020,000D)	1	Not allowed	From MWL
>Referenced Study Sequence	(0008,1110)	2	Not allowed	From MWL
>>Referenced SOP Class UID	(0008,1150)	1	Not allowed	From MWL

>>Referenced SOP Instance UID	(0008,1155)	1	Not allowed	From MWL
>Accession Number	(0008,0050)	2	Not allowed	From MWL
>Requested Procedure ID	(0040,1001)	2	Not allowed	From MWL
>Requested Procedure Code Sequence	(0032,1064)	3	Not allowed	From MWL
>Requested Procedure Description	(0032,1060)	2	Not allowed	From MWL
>Scheduled Procedure Step ID	(0040,0009)	2	Not allowed	From MWL
>Scheduled Procedure Step Description	(0040,0007)	2	Not allowed	From MWL / Entered by User
>Scheduled Protocol Code Sequence	(0040,0008)	2	Not allowed	From MWL / Entered by User
Patient's Name	(0010,0010)	2	Not allowed	From MWL / Entered by User
Patient ID	(0010,0020)	2	Not allowed	From MWL / Entered by User
Patient's Birth Date	0010,0030)	2	Not allowed	From MWL / Entered by User
Patient's Sex	(0010,0040)	2	Not allowed	From MWL / Entered by User
Referenced Patient Sequence	(0008,1120)	2	Not allowed	From MWL / Entered by User
Protocol Name	(0018,1030)	Not allowed	1	Updated at the end of Examination

Table 9.5-2: Performed Procedure Step Information Module attributes

Attribute Name	Tag	Req.Type N-CREATE	Req.Type N-SET	Note
Performed Procedure Step ID	(0040,0253)	1	Not allowed	Returned by MPPS server after Begin Study
Performed Station AE Title	(0040,0241)	1	Not allowed	Returned by MPPS server after Begin Study
Performed Station Name	(0040,0242)	2	Not allowed	Returned by MPPS server after Begin Study
Performed Location	(0040,0243)	2	Not allowed	Returned by MPPS server after Begin Study
Performed Procedure Step Start Date	(0040,0244)	1	Not allowed	Returned by MPPS server after Begin Study

Performed Procedure Step Start Time	(0040,0245)	1	Not allowed	Returned by MPPS server after Begin Study
Performed Procedure Step Status	(0040,0252)	1	3	This value can be: In Progress, Completed or Discontinued. Returned by MPPS server after Begin Study
Performed Procedure Step Description	(0040,0254)	2	3	Returned by MPPS server after Begin Study
Performed Procedures Type Description	(0040,0255)	2	3	Returned by MPPS server after Begin Study
Performed Procedure Code Sequence	(0008,1032)	2	3	Returned by MPPS server after Begin Study
>Code Value	(0008,0100)			
>Coding Scheme Designator	(0008,0102)			
>Code Meaning	(0008,0104)			
Performed Procedure Step End Date	(0040,0250)	3	2	Updated on End Study
Performed Procedure Step End Time	(0040,0251)	3	2	Updated on End Study

The following table specifies the attributes which describe the acquisition of Images during the Performance of the MPPS.

Table 9.5-3: Image Acquisition Results Module attributes

Attribute Name	Tag	Req.Type N-CREATE	Req.Type N-SET	Note
Modality	(0008,0060)	1	Not allowed	Type of Equipment =US. Retrieved from MWL / entered by user
Study ID	(0020,0010)	2	Not allowed	Retrieved from MWL / entered by user
Performed Protocol Code Sequence	(0040,0260)	2	2	Updated at the end study by MPPS
Performed Series Sequence	(0040,0340)	2	1	Updated at the end study by MPPS
>Performing Physician's Name	(0008,1050)	2	2	Updated at the end study by MPPS
>Protocol Name	(0018,1030)	1	1	Updated at the end study by MPPS

>Operators' Name	(0008,1070)	2	2	Updated at the end study by MPPS
>Series Instance UID	(0020,000E)	1	1	Updated at the end study by MPPS
>Series Description	(0008,103E)	2	2	Updated at the end study by MPPS
>Retrieve AE Title	(0008,0054)	2	2	Updated at the end study by MPPS
>Referenced Image Sequence	(0008,1140)	2	2	Updated at the end study by MPPS
>>Referenced SOP Class UID	(0008,1150)	1	1	Updated at the end study by MPPS
>>Referenced SOP Instance UID	(0008,1155)	1	1	Updated at the end study by MPPS
>Referenced Non Image composite SOP Instance Sequence	(0040,0220)	2	2	Updated at the end study by MPPS
>>Referenced SOP Class UID	(0008,1150)	1	1	Updated at the end study by MPPS
>>Referenced SOP Instance UID	(0008,1155)	1	1	Updated at the end study by MPPS

9.6 Overview of the Applied Storage Commitment IOD

This section specifies in detail the applied attributes in the N-ACTION service element of this supported SOP class.

Table 9.6-1: Storage Commitment Attribute Module

Attribute Name	Tag	Usage	Note
Transaction UID	(0008,1195)	M	Uniquely generated by the equipment
Retrieve AE Title	(0008,0054)	U	Not used
Storage Media File Set ID	(0088,0130)	U	Not used
Storage Media File Set UID	(0088,0140)	U	Not used
Referenced SOP Sequence	(0008,1199)	M	Supported
>Referenced SOP Class UID	(0008,1150)	M	Supported
>Referenced SOP Instance UID	(0008,1155)	M	Supported
>Retrieve AE Title	(0008,0054)	U	Not used
>Storage Media File Set	(0088,0130)	U	Not used

ID			
>Storage Media File – Set UID	(0088,0140)	U	Not used
Failed SOP Sequence	(0008,1198)	M	Supported N-Event-RQ
>References SOP Class UID	(0008,1150)	M	Supported
>Referenced SOP Instance UID	(0008,1155)	M	Supported
>Failure Reason	(0008,1197)	M	Supported

9.7 Overview of the Applied Comprehensive SR IOD

The Comprehensive SR IOD specifies a class of documents, the content of which may include textual and coded information, numeric measurement values, references to the DICOM Composite Instances.

The modules selected from the IOD module table of DICOM 3.0 are given in the table below.

Table 9.7-1: Applied Modules in the Comprehensive SR IOD

Information Entity	Module	Usage	Details
Patient	Patient	M	See Table 9.1-3
Study	General Study	M	See Table 9.1-4
	Patient Study	U	See Table 9.1.5
Series	SR Document Series	M	See Table 9.7-2
Equipment	General Equipment	M	See Table 9.1-7
Document	SR Document General	M	See Table 9.7-3
	SR Document Content	M	See Table 9.7-4
	SOP Common	M	See Table 9.7-5

The details of these applied Comprehensive SR IOD modules are given in the tables. The list of possible values is given, if applicable. If an attribute may be present conditionally/optionally or may contain a zero length value, this is also indicated.

Table 9.7-2: SR Document Series

Attribute Name	Tag	Type	Note
Modality	(0008,0060)	1	SR
SeriesInstanceUID	(0020,000E)	1	Generated at creation of the Series
SeriesNumber	(0020,0011)	1	Auto generated
Referenced Performed Procedure Step Sequence	(0008,1111)	2	Value set to EMPTY

Table 9.7-3: SR Document General

Attribute Name	Tag	Type	Note
Content Date	(0008,0023)	1	Date of SR document creation, based upon when user action that content creation started
Content Time	(0008,0033)	1	Time of SR document creation, based upon when user action that content creation started
Instance Number	(0020,0013)	1	See Table 9.1-8: General Image Module
Verifying Observer Sequence	(0040,A073)	1C	Not Used
Predecessor Documents Sequence	(0040,A360)	1C	Not Used
Referenced Request Sequence	(0040,A370)	1C	Not Used
Performed Procedure Code Sequence	(0040,A372)	2	Value set to EMPTY
Current Requested Procedure Evidence Sequence	(0040,A375)	1C	Not Used
Pertinent Other Evidence Sequence	(0040,A385)	1C	Not Used

Completion Flag	(0040,A491)	1	PARTIAL
Verification Flag	(0040,A493)	1	UNVERIFIED
Identical Documents Sequence	(0040,A525)	1C	Not Used

Table 9.7-4: SR Document Content

Attribute Name	Tag	Type	Note
Observation Date Time	(0040,A032)	1C	Not used
Value Type	(0040,A040)	1	CONTAINER
Concept Name Code Sequence	(0040,A043)	1C	Depending on the exam type - 'Vascular Ultrasound Procedure Report' 'Generic Ultrasound Structured Report'
Continuity Of Content	(0040,A050)	1C	SEPARATE
Content Template Sequence	(0040,A504)	1C	Depending on application template – Template ID's used are 'Vascular Ultrasound Procedure Reports' (TID 5100) and 'Generic Ultrasound Structured Report' (TID CS1000)
>Mapping Resource	(0008,0105)	1	Mapping Resource that defines the template. 1)'DCMR' for 'Vascular Ultrasound Procedure Report' 2)'CSCMR' for 'Generic Ultrasound Structured Report'
>Template Identifier	(0040,DB00)	1	Template identifier 1)'5100' for 'Vascular Ultrasound Procedure Report' 2)'CS1000' for 'Generic Ultrasound Structured Report'
Content Sequence	(0040,A730)	1C	Depending on application template - 'Vascular Ultrasound Procedure Report'

			'Generic Ultrasound Structured Report'
--	--	--	--

Table 9.7-5: SOP Common Module

Attribute Name	Tag	Type	Note
SOP Class UID	(0008,0016)	1	For Comprehensive SR: "1.2.840.10008.5.1.4.1.1.88.33"
SOP Instance UID	(0008,0018)	1	Generated when SR is created
Specific Character Set	(0008,0005)	1C	Set according to selected language on the scanner. See section 7.

9.7.1 Overview of SR Document Content Descriptions

9.7.1.1 SR Document Content Template

Carestream Ultrasound Systems supports the following DICOM SR root templates for Comprehensive SR SOP Instances created, generated, or displayed by the systems.

Table 9.7-6: DICOM SR Root Templates

Template ID	Template Name
5100	Vascular Ultrasound Procedure Reports
CS1000	Generic Ultrasound Structured Report

9.7.1.2 DICOM Standard Extended and Private Template

Carestream Ultrasound Systems supports the DICOM standard extended and private defined templates in the following sections.

9.7.1.2.1 TID 5100 Vascular SR - DICOM Standard Extended

DICOM Standard Extended Table 9.7-7 TID 5100 Vascular Ultrasound Procedure Reports for user selected Exam Type Cerebro Vascular or Peripheral Vascular for Carestream Ultrasound Systems.

Table 9.7-7: TID 5100 - VASCULAR ULTRASOUND REPORT

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (125100, DCM, "Vascular Ultrasound Procedure Report")	1	M		
2	>	HAS OBS CONTEXT	INCLUDE	DTID (1001) Observation Context	1	M		Person Name value of Performing-Physician
3	>	CONTAINS	CONTAINER	EV (111028, DCM, "Image Library")	1	U		
4	>>	CONTAINS	IMAGE	No purpose of reference	1-n	M		DCID(12224) Ultrasound Image Modes
5	>	CONTAINS	INCLUDE	DTID (5102) Vascular Procedure Summary Section	1	U		
6	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV(G-A101, SRT, "Left") \$Anatomy = DCID (12105) Intracranial Cerebral Vessels
7	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12105) Intracranial Cerebral Vessels
8	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-40501, SRT, "Blood Vessel of Head") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12106) Intracranial Cerebral Vessels (unilateral)
9	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U		\$SectionScope = DT (T-45005, SRT, "Artery of neck") \$SectionLaterality = EV(G-A101, SRT, "Left") \$Anatomy = DCID (12104) Extracranial Arteries \$AnatomyRatio = DCID

							(12123) Carotid Ratios
10	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-45005, SRT, "Artery of neck") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12104) Extracranial Arteries \$AnatomyRatio = DCID (12123) Carotid Ratios
11	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-47040, SRT, "Artery of Lower Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12109) Lower Extremity Arteries
12	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-47040, SRT, "Artery of Lower Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12109) Lower Extremity Arteries
13	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-47020, SRT, "Artery Of Upper Extremity") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12107) Upper Extremity Arteries
14	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-47020, SRT, "Artery

							Of Upper Extremity") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12107) Upper Extremity Arteries
15	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") \$SectionLaterality = EV (G-A101, SRT, "Left") \$Anatomy = DCID (12115) Renal Vessels \$AnatomyRatio = DCID (12124) Renal Ratios
16	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-71019, SRT, "Vascular Structure Of Kidney") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12115) Renal Vessels \$AnatomyRatio = DCID (12124) Renal Ratios
17	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-46002, SRT, "Artery of Abdomen") A101, SRT, "Left") \$Anatomy = DCID (12111) Abdominal Arteries (lateral)
18	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound	1	U	\$SectionScope = DT (T-

				Section			46002, SRT, "Artery of Abdomen") \$SectionLaterality = EV (G-A100, SRT, "Right") \$Anatomy = DCID (12111) Abdominal Arteries (lateral)
19	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-46002, SRT, "Artery of Abdomen") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12112) Abdominal Arteries (unilateral)
20	>	CONTAINS	INCLUDE	DTID (5103) Vascular Ultrasound Section	1	U	\$SectionScope = DT (T-487A0, SRT, "Vein of Abdomen") \$SectionLaterality = EV (G-A103, SRT, "Unilateral") \$Anatomy = DCID (12114) Abdominal Veins (unilateral)

Table 9.7-8: TID 5102 - Vascular Procedure Summary Section

Description:								
Comments and observations of the procedure of immediate clinical interest								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (121111, DCM, "Summary")	1	M		
2	>	CONTAINS	TEXT	DCID (12101) Vascular Summary	1-n	M		Text Value of Report Remarks

Table 9.7-9: TID 5103 - Vascular Ultrasound Section

Description:								
Sections of a vascular ultrasound report are section containers of an anatomic region consisting of measurement group containers that contain the measurements								
Parameters:								
<p>\$SectionScope: The concept name of the section heading modifier</p> <p>\$SectionLaterality: The laterality (if any) of the anatomy in this section heading</p> <p>\$Anatomy: The concept name of the vascular anatomy</p> <p>\$AnatomyRatio: The concept name of anatomy-coordinated ratio concepts</p>								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (121070, DCM, "Findings")	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C0E3, SRT, "Finding Site")	1	M		\$SectionScope
3	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, "Laterality")	1	M		\$SectionLaterality
4	>	CONTAINS	INCLUDE	DTID (5104) Vascular Measurement Group	1-n	M		\$AnatomyGroup = \$Anatomy
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement = \$AnatomyRatio

Table 9.7-10: TID 5104 - Vascular Measurement Group

Description:								
This template is an anatomy specific container of measurements								
Parameters:								
\$AnatomyGroup: The concept name of the vascular anatomy								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$AnatomyGroup	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-A1F8, SRT, "Topographical Modifier")	1	U		DCID (12116) Vessel Segment Modifiers
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	M		\$Measurement = DCID (12119) Vascular Ultrasound Property

Table 9.7-11: CID 12104 - Extracranial Arteries

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-45170	Carotid Bulb
SRT	T-45100	Common Carotid Artery
SRT	T-45200	External Carotid Artery
SRT	T-45300	Internal Carotid Artery
SRT	T-46100	Subclavian Artery
SRT	T-45700	Vertebral Artery

Table 9.7-12: CID 12105 - Intracranial Arteries

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-45540	Anterior Cerebral Artery
SRT	T-45308	Carotid Siphon
SRT	T-45300	Internal Carotid Artery
SRT	R-102BD	Terminal internal carotid artery
SRT	R-1024F	Middle Cerebral Artery M1 Segment

SRT	R-10251	Middle Cerebral Artery M2 Segment
SRT	T-45400	Ophthalmic Artery
SRT	R-10253	Posterior Cerebral Artery P1 Segment
SRT	R-10255	Posterior Cerebral Artery P2 Segment

Table 9.7-13: CID 12106 - Intracranial Arteries (unilateral)

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-45800	Basilar Artery

Table 9.7-14: CID 12107 – Upper Extremity Arteries

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-47100	Acillary Artery
SRT	T-47160	Brachial Artery
SRT	T-46010	Innominate Artery
SRT	T-47300	Radial Artery
SRT	T-46100	Subclavian Artery
SRT	T-47200	Ulnar Artery

Table 9.7-15: CID 12109 - Lower Extremity Arteries

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-47440	Profunda Femoris Artery
SRT	T-47400	Common Femoral Artery
SRT	T-47500	Popliteal Artery
SRT	T-46710	Common Iliac Artery
SRT	T-46910	External Iliac Artery
SRT	T-47600	Posterior Tibial Artery
SRT	T-47741	Dorsalis Pedis Artery
SRT	T-47403	Superficial Femoral Artery
SRT	T-47700	Anterior Tibial Artery

SRT	T-47630	Peroneal Artery
-----	---------	-----------------

Table 9.7-16: CID 12111 - Abdominal Arteries (lateral)

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-46980	Ovarian Artery
SRT	T-46820	Uterine Artery
SRT	T-46710	Common Iliac Artery

Note: Adding SRT, T-46710, Common Iliac Artery of this template is a private extension

Table 9.7-17: CID 12112 - Abdominal Arteries (unilateral)

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-42000	Aorta
SRT	T-46400	Celiac Axis
SRT	T-46421	Common Hepatic Artery
SRT	T-46422	Proper Hepatic Artery
SRT	T-46460	Splenic Artery
SRT	T-46510	Superior Mesenteric Artery
SRT	T-F1810	Umbilical Artery

Note: Adding SRT, TF1810, Umbilical Artery under CID 12112 instead of CID 12111 is a private extension. No attempt is made to differentiate between the two umbilical arteries.

Table 9.7-18: CID 12114 - Abdominal Veins (unilateral)

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-48710	Inferior Vena Cava
SRT	T-48727	Left Hepatic Vein
SRT	T-48726	Middle Hepatic Vein
SRT	T-48725	Right Hepatic Vein
SRT	T-48810	Portal Vein
SRT	T-4881F	Left Main Branch of Portal Vein

SRT	T-4882A	Right Main Branch of Portal Vein
SRT	T-48890	Splenic Vein
SRT	T-48840	Superior Mesenteric Vein

Table 9.7-19: CID 12115 – Renal Vessels

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-46600	Renal Artery

Table 9.7-20: CID 12116 - Vessel Segment Modifiers

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	G-A119	Distal
SRT	G-A188	Mid-longitudinal
SRT	G-A118	Proximal

Table 9.7-21: CID 12119 - Vascular Ultrasound Property

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
INCLUDE CID 12120 Blood Velocity Measurements		
INCLUDE CID 12121 Vascular Indices and Ratios		
INCLUDE CID 12122 Other Vascular Properties		

Table 9.7-22: CID 12120 - Blood Velocity Measurements

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11653-3	End Diastolic Velocity
LN	11726-7	Peak Systolic Velocity
LN	11665-7	Minimum Diastolic Velocity
LN	20351-1	Time averaged mean velocity

LN	11692-1	Time averaged peak velocity
----	---------	-----------------------------

Table 9.7-23: CID 12121 – Vascular Indices and Ratios

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	12008-9	Pulsatility Index
LN	12023-8	Resistivity Index
LN	12144-2	Systolic to Diastolic Velocity Ratio
LN	20351-1	Time averaged mean velocity
LN	11692-1	Time averaged peak velocity

Table 9.7-24: CID12122 - Other Vascular Properties

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	20168-1	Acceleration Time
LN	20217-6	Deceleration Time

Table 9.7-25: CID12123 - Carotid Ratios

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	33868-1	ICA/CCA velocity ratio

Table 9.7-26: CID12124 – Renal Ratios

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	33869-9	Renal Artery/Aorta velocity ratio

Table 9.7-27: CID 12224 - Ultrasound Image Modes

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	G-03A2	2D mode

Table 9.7-28: TID 300 - Measurement

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			NUM	\$Measurement	1	M		Refer CID 82 - Units of Measurement. Note: Peak Systole uses cm/s.

9.7.1.2.2 TID CS1000 Generic Ultrasound Structured Report – Private User Defined

The Generic Ultrasound Structured Report template Table 9.7-28 provides the specific codes for user defined measurements and calculations with public/private code designator schemes.

The TID CS1000 is used when any ExamType other than Cerebro Vascular and Peripheral Vascular is active.

Table 9.7-29: TID CS1000 – Generic Ultrasound Structured Report

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	EV (PC-00001, PC, "Generic Ultrasound Structured Report")	1	M		
2	>	HAS OBS CONTEXT	INCLUDE	DTID (1001) Observation Context	1	M		Person Name value of Performing-Physician
3	>	CONTAINS	CONTAINER	EV (111028, DCM, "Image Library")	1	U		
4	>>	CONTAINS	IMAGE	No purpose of reference	1-n	M		DCID(12224) Ultrasound Image Modes
5	>	CONTAINS	INCLUDE	DTID (CS100) Summary Section	1	U		
6	>	CONTAINS	INCLUDE	DTID (CS200) Ultrasound Measurements Section	1	U		

Table 9.7-30: TID CS100 - Summary Section

Description:								
Comments and observations of the procedure of immediate clinical interest								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (121111, DCM, "Summary")	1	M		
2	>	CONTAINS	TEXT	EV(DCM,121106,"Comment")	1	M		Text Value of Report Remarks

Table 9.7-31: TID CS200 – Ultrasound Measurements Section

Description:								
Main section of the Generic Ultrasound Structured Report. Consists of a number of containers relating to specific examination types.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00002, PC, “Ultrasound Measurements”)	1	U		
2	>	CONTAINS	INCLUDE	DTID (CS1300) Abdominal Measurements	1	U		
3	>	CONTAINS	INCLUDE	DTID (CS1400) Breast Measurements	1-2	U		\$Laterality = DCID (CS3102) Side (Breast Identification)
4	>	CONTAINS	INCLUDE	DTID (CS1500) Early Ob Measurements	1	U		
5	>	CONTAINS	INCLUDE	DTID (CS1510) Ob Measurements	1	U		
6	>	CONTAINS	INCLUDE	DTID (CS1520) Ob Fetal Measurements	1	U		
7	>	CONTAINS	INCLUDE	DTID (CS1600) Gyn Measurements	1	U		
8	>	CONTAINS	INCLUDE	DTID (CS1700) Prostate Measurements	1	U		
9	>	CONTAINS	INCLUDE	DTID (CS1800) Renal Measurements	1	U		
10	>	CONTAINS	INCLUDE	DTID (CS1900) Testis Measurements	1	U		
11	>	CONTAINS	INCLUDE	DTID (CS2000) Thyroid Measurements	1	U		

Table 9.7-32: TID CS1300 – Abdominal Measurements

Description:								
Container with measurements related to the Abdominal Exam Type.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00011, PC,				

				“Abdominal Measurements”)				
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	M		\$Measurement = DCID (CS3100) Abdominal Ultrasound Property
3	>	CONTAINS	INCLUDE	DTID (CS3011) Kidney Measurement Group	1	U		
4	>	CONTAINS	INCLUDE	DTID (CS3000) Volume Measurement from LHW	1	U		\$Anatomy = DT (TC3000, SRT, “Spleen”)
5	>	CONTAINS	INCLUDE	DTID (CS3001) Vessel Diameter Measurement Group	1	U		\$Vessel = DT (T-42000, SRT, “Aorta”)
6	>	CONTAINS	INCLUDE	DTID (CS3001) Vessel Diameter Measurement Group	1	U		\$Vessel = DT (T-46710, SRT, “Common Iliac Artery”)
7	>	CONTAINS	CONTAINER	DT (T-48710, SRT, “Inferior Vena Cava”	1	U		
8	>>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (18006-7, LN, “Inferior Vena Cava Diameter)
9	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1-n	U		\$MeasurmentGroup = DCID (CS3101) Abdominal Vessels for Generic Report

Table 9.7-33: TID CS1400 – Breast Measurements

Description:								
Container with measurements related to the Breast Exam Type. Up to two containers can be included. Each container has measurements related to one breast only.								
Parameters:								
\$Laterality: Identifies the breast for which this container has measurements								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00012, PC, “Breast Measurements”)				
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	U		\$Laterality
3	>	CONTAINS	INCLUDE	DTID (CS3003) Breast Mass Measurement Group	1-n	U		\$Mass = DCID (CD3103) Mass Identification

Table 9.7-34: TID CS1500 – Early Ob Measurements

Description:								
Container with measurements related to the Early Ob Exam Type (fetus properties are not here, but in TID CS1520).								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00013, PC, “Early Ob Measurements”)				
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement=DCID (CS3104) Early Ob Property
3	>	CONTAINS	INCLUDE	DTID (CS3004) Ovary Measurement Group	1-2	U		\$Laterality = “Left” or “Right”
4	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1-n	U		\$AnatomyGroup = DCID (CS3105) Arteries for Early Ob

Table 9.7-35: TID CS1510 – Ob Measurements

Description:								
Container with measurements related to the Ob Exam Type (fetus properties are not here, but in TID CS1520).								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00014, PC, “Ob Measurements”)				
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement=DCID (CS3106) Ob Property
3	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1-n	U		\$AnatomyGroup = DCID (CS3107) Arteries for Ob

Table 9.7-36: TID CS1520 – Ob Fetal Measurements

Description:								
Container with fetus property measurements related to the Ob and Early Ob Exam Types								

	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00015, PC, “Ob Fetal Measurements”)				
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Measurement=DCID (CS3108) Ob Fetal Conclusions
3	>	CONTAINS	INCLUDE	DTID (CS3005) Fetal Biometry Measurement Group	1-n	U		\$Biometry = DCID (CS3109) Ob Fetal Biometry for Generic Report
4	>	CONTAINS	INCLUDE	DTID (CS3006) Fetal Biometry Ratio Calculation Group	1-n	U		\$Biometry = DCID (12004) Fetal Biometry Ratios
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement=DT (8867-4, LN, “Heart Rate”)
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-4	U		\$Measurement=DCID (12008) OB-GYN Amniotic Sac
7	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1	U		\$AnatomyGroup = DT (T-F1810, SRT, “Umbilical Artery”)

Table 9.7-37: TID CS1600 – Gyn Measurements

Description:								
Container with measurements related to the Gyn Exam Type								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00016, PC, “Gyn Measurements”)				
2	>	CONTAINS	INCLUDE	DTID (CS3008) Uterus Measurement Group	1	U		
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	U		\$Biometry = DCID (CS3110) Gyn Ultrasound Property
4	>	CONTAINS	INCLUDE	DTID (CS3004) Ovary Measurement Group	1-2	U		\$Laterality = “Left” or “Right”
5	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder Volume Measurement from LHW	1	U		\$BladderCondition=DT (PC-13000, PC, “Pre Void Bladder”)
6	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder	1	U		\$BladderCondition=DT (PC-

				Volume Measurement from LHW				13001, PC, "Post Void Bladder")
7	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement=DT (PC-13002, PC, "Micturated Volume")
7	>	CONTAINS	INCLUDE	DTID (CS3009) Follicle Volume Measurement from Diameter	1	U		\$Follicle=DCID (CS3111) Follicle Identification

Table 9.7-38: TID CS1700 – Prostate Measurements

Description:								
Container with measurements related to the Prostate Exam Type								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00017, PC, "Prostate Measurements")				
2	>	CONTAINS	INCLUDE	DTID (CS3010) Prostate Volume Measurement from LHW	1	U		
3	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder Volume Measurement from LHW	1	U		\$BladderCondition=DT (PC-13000, PC, "Pre Void Bladder")
4	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder Volume Measurement from LHW	1	U		\$BladderCondition=DT (PC-13001, PC, "Post Void Bladder")
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement=DT (PC-13002, PC, "Micturated Volume")

Table 9.7-39: TID CS1800 – Renal Measurements

Description:								
Container with measurements related to the Renal Exam Type								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00018, PC, "Renal Measurements")				

2	>	CONTAINS	INCLUDE	DTID (CS3011) Kidney Measurement Group	1	U		
3	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder Volume Measurement from LHW	1	U		\$BladderCondition=DT (PC-13000, PC, "Pre Void Bladder")
4	>	CONTAINS	INCLUDE	DTID (CS3008) Bladder Volume Measurement from LHW	1	U		\$BladderCondition=DT (PC-13001, PC, "Post Void Bladder")
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement=DT (PC-13002, PC, "Micturated Volume")
6	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1-n	U		\$MeasurementGroup = DT (T-46600, SRT, "Renal Artery")
7	>	CONTAINS	INCLUDE	DTID (CS3002) Vascular Measurement Group for Generic Report	1-n	U		\$MeasurementGroup = DT (T-42000, SRT, "Aorta")
8	>	CONTAINS	INCLUDE	DTID (CS3012) Renal Ratio	1	U		\$Laterality = EV (G-A100, SRT, "Right")
9	>	CONTAINS	INCLUDE	DTID (CS3012) Renal Ratio	1	U		\$Laterality = EV (G-A101, SRT, "Left")

Table 9.7-40: TID CS1900 – Testis Measurements

Description:								
Container with measurements related to the Testis Exam Type								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00019, PC, "Testicular Measurements")				
2	>	CONTAINS	INCLUDE	DTID (CS3013) Volume Measurement from LHW (lateral)	1	U		\$Anatomy=DT (T-94000, SRT, "Testis") \$Laterality = EV (G-A100, SRT, "Right")
3	>	CONTAINS	INCLUDE	DTID (CS3013) Volume Measurement from LHW (lateral)	1	U		\$Anatomy=DT (PC-14000, SRT, "Epididymis") \$Laterality = EV (G-A100, SRT, "Right")
4	>	CONTAINS	INCLUDE	DTID (CS3014) Scrotum Measurement Group	1	U		\$Laterality = EV (G-A100, SRT, "Right")

5	>	CONTAINS	INCLUDE	DTID (CS3013) Volume Measurement from LHW (lateral)	1	U		\$Anatomy=DT (T-94000, SRT, "Testis") \$Laterality = EV (G-A101, SRT, "Left")
6	>	CONTAINS	INCLUDE	DTID (CS3013) Volume Measurement from LHW (lateral)	1	U		\$Anatomy=DT (PC-14000, SRT, "Epididymis") \$Laterality = EV (G-A101, SRT, "Left")
7	>	CONTAINS	INCLUDE	DTID (CS3014) Scrotum Measurement Group	1	U		\$Laterality = EV (G-A101, SRT, "Left")

Table 9.7-41: TID CS2000 – Thyroid Measurements

Description:								
Container with measurements related to the Thyroid Exam Type								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00020, PC, "Thyroid Measurements")				
2	>	CONTAINS	INCLUDE	DTID (CS3015) Thyroid Measurement Group	1-n	U		\$ThyroidAnatomy=DCID (CS3112) Thyroid Feature \$Laterality = EV (G-A100, SRT, "Right")
3	>	CONTAINS	INCLUDE	DTID (CS3015) Thyroid Measurement Group	1-n	U		\$ThyroidAnatomy=DCID (CS3112) Thyroid Feature \$Laterality = EV (G-A101, SRT, "Left")
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement=DT (PC-15000, PC, "Thyroid Isthmus")

Table 9.7-42: TID CS3000 – Volume Measurement from LHW

Description:								
Container which groups an LHW Volume measurement for a given anatomy.								
Parameters:								
\$Anatomy: The anatomy being measured.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Anatomy	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A22A, SRT, "Length")
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (121207, DCM, "Height")
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A220, SRT, "Width")
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 2-4 are all present	\$Measurement = DT (G-D705, SRT, "Volume")

Table 9.7-43:TID CS3001- Vessel Diameter Measurement Group

Description:								
Container which groups Anterior-Posterior and Transverse Diameter measurements for a vessel.								
Parameters:								
\$Anatomy: The vessel being measured.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Vessel	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	U		DCID (244) Laterality
3	>	HAS CONCEPT MOD	CODE	EV (G-A1F8, SRT, “Topographical Modifier”)	1	U		DCID (12116) Vessel Segment Modifiers
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (PC-00500, PC, “Anterior-Posterior Diameter”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (PC, PC-00501, “Transverse Diameter”)

Table 9.7-44: TID CS3002 - Vascular Measurement Group for Generic Report

Description:								
Container which groups vascular measurements for a given abdominal vascular anatomy.								
Parameters:								
\$AnatomyGroup: The concept name of the vascular anatomy								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$AnatomyGroup	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	U		DCID (244) Laterality
3	>	HAS CONCEPT MOD	CODE	EV (G-A1F8, SRT, “Topographical Modifier”)	1	U		DCID (12116) Vessel Segment Modifiers
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1-n	M		\$Measurement = DCID (12119) Vascular Ultrasound Property

Table 9.7-45: TID CS3003 – Breast Mass Measurement Group

Description:								
Container which groups an LHW Volume measurement for a given breast mass.								
Parameters:								
\$Mass: Identifier for the mass being measured.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Anatomy	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (PC-12010, PC, “Mass Id”)	1	M		\$Mass
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A22A, SRT, “Length”)
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (121207, DCM, “Height”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A220, SRT, “Width”)
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 2-5 are all present	\$Measurement = DT (G-D705, SRT, “Volume”)

Table 9.7-46: TID CS3004 – Ovary Measurement Group

Description:								
Container which groups an LHW Volume measurement for a given ovary								
Parameters:								
\$Laterality: Identifies the ovary								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (T-87000, SRT, “Ovary”)	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF \$Laterality = “Right” and Rows 3-5 are all present	\$Measurement = DT (12165-7, LN, “Right Ovary Volume”)
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Right”	\$Measurement = DT (11841-4, LN, “Right Ovary Length”)
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Right”	\$Measurement = DT (11858-8, LN, “Right Ovary Height”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Right”	\$Measurement = DT (11830-7, LN, “Right Ovary Width”)
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF \$Laterality = “Left” and Rows 7-9 are all present	\$Measurement = DT (12164-0, LN, “Left Ovary Volume”)
7	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Left”	\$Measurement = DT (11840-6, LN, “Left Ovary Length”)
8	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Left”	\$Measurement = DT (11857-0, LN, “Left Ovary Height”)
9	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U	IFF \$Laterality = “Left”	\$Measurement = DT (11829-9, LN, “Left Ovary Width”)

Table 9.7-47: TID CS3005 – Fetal Biometry Measurement Group

Description: Container which groups a fetal biometry measurement with calculations of gestational age and growth percentile rank (if these calculations are available)								
Parameters: \$Biometry: A fetal biometry								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (125005, DCM, “Biometry Group”)	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = \$Biometry
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (18185-9, LN, “Gestational Age”)
4	>	CONTAINS	TEXT	DT (PC-20000, PC, “Gestational Age Author”)	1	MC	IFF Row 3 is present	Name of the author of the calculation used for gestational age
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (125012, DCM, “Growth Percentile Rank”)
6	>	CONTAINS	TEXT	DT (PC-20001, PC, “Growth Percentile Rank Author”)	1	MC	IFF Row 5 is present	Name of the author of the calculation used for growth percentile rank.

Table 9.7-48: TID CS3006 – Fetal Biometry Ratio Calculation Group

Description: Container which groups calculations of a fetal biometry ratio with calculations of the min and max percentiles based on the average ultrasound age (if available). Parameters: \$BiometryRatio: A fetal biometry ratio								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (125005, DCM, “Biometry Group”)	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = \$BiometryRatio
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (PC-20002, PC, “Ratio min percentile based on Average US Age”)
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (PC-20003, PC, “Ratio max percentile based on Average US Age”)
5	>	CONTAINS	TEXT	DT (PC-20004, PC, “Ratio percentile Author”)	1	MC	IFF Row 3 and/or Row 4 is present	Name of the author of the calculation used for the min and max percentiles.

Table 9.7-49: TID CS3007 – Uterus Measurement Group

Description:								
Container which groups an LHW volume measurement for the uterus								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (T-83000, SRT, "Uterus")	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (11842-2, LN, "Uterus Length")
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (11859-6, LN, "Uterus Height")
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (11865-3, LN, "Uterus Width")
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M	IFF Rows 2-4 are all present	\$Measurement = DT (33192-6, LN, "Uterus Volume")

Table 9.7-50: TID CS3008 – Bladder Volume Measurement from LHW

Description:								
Container which groups an LHW volume measurement for the pre or post void bladder								
Parameters:								
\$BladderCondition: Pre or Post void.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$BladderCondition	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15300-7, LN, “Urinary bladder Length by US”)
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15298-3, LN, “Urinary bladder Height by US”)
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15299-1, LN, “Urinary bladder Width by US”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 2-4 are all present	\$Measurement = DT (15309-8, LN, “Urinary bladder Volume from HWL by US”)

Table 9.7-51: TID CS3009 – Follicle Volume Measurement from Diameter

Description:								
Container which groups a diameter-based volume measurement for a follicle								
Parameters:								
\$Follicle: Identifies the follicle.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Follicle	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (11793-7, LN, "Follicle Diameter")
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (G-D705, SRT, "Volume")

Table 9.7-52: TID CS3010 – Prostate Volume Measurement from LHW

Description:								
Container which groups an LHW volume measurement for the prostate								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (T-9200B, SRT, Prostate)	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15303-1, LN, "Prostate Length by US")
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15301-5, LN, "Prostate Height by US")
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15302-3, LN, "Prostate Width by US")

5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 2-4 are all present	\$Measurement = DT (15308-0, LN, "Prostate Volume from HWL by US")
---	---	----------	---------	------------------------	---	----	------------------------------	--

Table 9.7-53: TID CS3011 – Kidney Measurement Group

Description:								
Container which groups kidney measurements (for both right and left kidney)								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (T-9200B, SRT, Prostate)	1	M		
2	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15291-8, LN, "Right Kidney Length")
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15289-2, LN, "Right Kidney Height")
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15290-0, LN, "Right Kidney Width")
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 2-4 are all present	\$Measurement = DT (15305-6, LN, "Right Kidney Volume")
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15288-4, LN, "Left Kidney Length")
7	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15286-8, LN, "Left Kidney Height")
8	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (15287-6, LN, "Left Kidney Width")
9	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 6-8 are all	\$Measurement = DT (15304-9, LN, "Left Kidney

							present	Volume”)
--	--	--	--	--	--	--	---------	----------

Table 9.7-54: TID CS3012 – Renal Ratio

<p>Description:</p> <p style="text-align: center;">Container which contains the calculation result of a renal ratio.</p> <p>Parameters:</p> <p style="text-align: center;">\$Laterality: Laterality of the kidney for which the renal ratio is calculated</p>								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$AnatomyGroup	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	M		\$Laterality
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (33869-9, LN, “Renal Artery/Aorta velocity ratio”)

Table 9.7-55: TID CS3013 – Volume Measurement from LHW (lateral)

Description: Container which groups an LHW Volume measurement for a given anatomy (lateral). Parameters: \$Anatomy: The anatomy being measured. \$Laterality: The laterality of the anatomy.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	\$Anatomy	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	M		\$Laterality
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A22A, SRT, “Length”)
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (121207, DCM, “Height”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A220, SRT, “Width”)
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 3-5 are all present	\$Measurement = DT (G-D705, SRT, “Volume”)

Table 9.7-56: TID CS3014 – Scrotum Measurement Group

Description: <p style="text-align: center;">Container which contains the measurements of one side of the scrotum.</p> Parameters: <p style="text-align: center;">\$Laterality: Laterality of scrotum measurements</p>								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (T-98000, SRT, “Scrotum”)	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	M		\$Laterality
3	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	M		\$Measurement = DT (PC-14001, PC, “Skin Thickness”)

Table 9.7-57: TID CS3015 – Thyroid Measurement Group

Description: Container which groups an LHW Volume measurement for a thyroid feature.								
Parameters: \$ThyroidAnatomy: The thyroid anatomy being measured. \$Laterality: The laterality of the anatomy.								
	NL	Rel with Parent	VT	Concept Name	VM	Req Type	Condition	Value Set Constraint
1			CONTAINER	DT (PC-00020, PC, “Thyroid Measurements”)	1	M		
2	>	HAS CONCEPT MOD	CODE	EV (G-C171, SRT, “Laterality”)	1	M		\$Laterality
3	>	HAS CONCEPT MOD	CODE	EV (G-A1F8, SRT, “Topographical Modifier”)	1	M		\$ThyroidAnatomy
4	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A22A, SRT, “Length”)
5	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (121207, DCM, “Height”)
6	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	U		\$Measurement = DT (G-A220, SRT, “Width”)
7	>	CONTAINS	INCLUDE	DTID (300) Measurement	1	MC	IFF Rows 4-6 are all present	\$Measurement = DT (G-D705, SRT, “Volume”)

Table 9.7-58: CID 244 - Laterality

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	G-A100	Right
SRT	G-A101	Left

SRT	G-A103	Unilateral
-----	--------	------------

Table 9.7-59: CID CS3100 – Abdominal Ultrasound Property

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PC	PC-11000	Liver Length
PC	PC-11001	Gall Bladder Wall
PC	PC-11002	Common Duct

Table 9.7-60: CID CS3101 – Abdominal Vessels for Generic Report

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-42000	Aorta
SRT	T-46710	Common Iliac Artery
SRT	T-46400	Celiac Axis
SRT	T-46421	Common Hepatic Artery
SRT	T-46460	Splenic Artery
SRT	T-46510	Superior Mesenteric Artery
SRT	T-48710	Inferior Vena Cava
SRT	T-48727	Left Hepatic Vein
SRT	T-48726	Middle Hepatic Vein
SRT	T-48725	Right Hepatic Vein
SRT	T-48810	Portal Vein
SRT	T-4881F	Left Main Branch of Portal Vein
SRT	T-4882A	Right Main Branch of Portal Vein
SRT	T-48890	Splenic Vein
SRT	T-48840	Superior Mesenteric Vein

Table 9.7-61: CID CS3102 – Side (Breast Identification)

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-04030	Left breast
SRT	T-04020	Right breast

Table 9.7-62: CID CS3103 – Mass Identification

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PC	PC-12001	Mass 1
PC	PC-12002	Mass 2
PC	PC-12003	Mass 3
PC	PC-12004	Mass 4
PC	PC-12005	Mass 5
PC	PC-12006	Mass 6
PC	PC-12007	Mass 7
PC	PC-12008	Mass 8

Table 9.7-63: CID CS3104 – Early Ob Property

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11961-0	Cervix Length
PC	PC-20005	Right Corpus Luteum
PC	PC-20006	Left Corpus Luteum

Table 9.7-64: CID CS3105 – Arteries for Early Ob

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-46980	Ovarian Artery
SRT	T-46820	Uterine Artery

Table 9.7-65: CID CS3106 – Ob Property

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11961-0	Cervix Length

Table 9.7-66: CID CS3107 – Arteries for Ob

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SRT	T-46980	Ovarian Artery
SRT	T-46820	Uterine Artery

Table 9.7-67: CID CS3108 – Ob Fetal Conclusions

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11885-1	Gestational Age by LMP
LN	11779-6	EDD from LMP
LN	11955-2	LMP
LN	11884-4	Average Ultrasound Age
LN	11781-2	EDD from average ultrasound age
LN	11727-5	Estimated Weight
PC	PC-20007	Estimated Fetal Weight Author
LN	11767-1	EFW percentile rank
PC	PC-20008	EFW Percentile Rank Author

Table 9.7-68: CID CS3109 – Ob Fetal Biometry for Generic Report

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
INCLUDE CID 12005 Fetal Biometry Measurements		
INCLUDE CID 12006 Fetal Long Bones Measurements		
INCLUDE CID 12007 Fetal Cranium		
INCLUDE CID 12009 Early Gestration Biometry Measurements		

Table 9.7-69: CID 12004 Fetal Biometry Ratios

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11947-9	HC/AC

LN	11871-1	FL/AC
LN	11872-9	FL/BPD
LN	11823-2	Cephalic Index
LN	11873-7	FL/HC

Table 9.7-70: CID 12005 Fetal Biometry Measurements

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11979-2	Abdominal Circumference
LN	11818-2	Anterior-Posterior Abdominal Diameter
LN	11820-8	Biparietal Diameter
LN	11963-6	Femur Length
LN	11965-1	Foot Length
LN	11984-2	Head Circumference
LN	11851-3	Occipital-Frontal Diameter
LN	11862-0	Transverse Abdominal Diameter
LN	11863-8	Trans Cerebellar Diameter

Table 9.7-71: CID 12006 Fetal Long Bones Measurements

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11966-2	Humerus length
LN	11967-7	Radius length
LN	11969-3	Ulna length
LN	11968-5	Tibia length
LN	11964-4	Fibula length
LN	11962-8	Clavicle length

Table 9.7-72: CID 12007 Fetal Cranium

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	12717-5	Lateral Ventricular Width
LN	11860-4	Cisterna Magna length
LN	12146-7	Nuchal Fold Thickness

LN	11629-3	Outer Orbital Diameter
LN	11863-8	Trans Cerebellar Diameter

Table 9.7-73: CID 12008 OB-GYN Amniotic Sac

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11624-4	First Quadrant Diameter
LN	11626-9	Second Quadrant Diameter
LN	11625-1	Third Quadrant Diameter
LN	11623-6	Fourth Quadrant Diameter
LN	11627-7	Amniotic Fluid Index

Note: Adding LN, 11627-7, Amniotic Fluid Index is a private extension

Table 9.7-74: CID 12009 Early Gestation Biometry Measurements

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	11957-8	Crown Rump Length
LN	11850-5	Gestational Sac Diameter
LN	11816-6	Yolk Sac length
LN	33069-6	Nuchal Translucency

Table 9.7-75: CID CS3110 – Gyn Ultrasound Property

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	12145-9	Endometrium Thickness
LN	11961-0	Cervix Length

Table 9.7-76: CID CS3111 – Follicle Identification

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PC	PC-13011	Right Ovarian Follicle 1
PC	PC-13012	Right Ovarian Follicle 2
PC	PC-13013	Right Ovarian Follicle 3

PC	PC-13014	Right Ovarian Follicle 4
PC	PC-13015	Right Ovarian Follicle 5
PC	PC-13016	Right Ovarian Follicle 6
PC	PC-13017	Right Ovarian Follicle 7
PC	PC-13018	Right Ovarian Follicle 8
PC	PC-13019	Right Ovarian Follicle 9
PC	PC-13020	Right Ovarian Follicle 10
PC	PC-13021	Left Ovarian Follicle 1
PC	PC-13022	Left Ovarian Follicle 2
PC	PC-13023	Left Ovarian Follicle 3
PC	PC-13024	Left Ovarian Follicle 4
PC	PC-13025	Left Ovarian Follicle 5
PC	PC-13026	Left Ovarian Follicle 6
PC	PC-13027	Left Ovarian Follicle 7
PC	PC-13028	Left Ovarian Follicle 8
PC	PC-13029	Left Ovarian Follicle 9
PC	PC-13030	Left Ovarian Follicle 10

Table 9.7-77: CID CS3112 – Thyroid Feature

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
PC	PC-15010	Lobe
PC	PC-15011	Nodule 1
PC	PC-15012	Nodule 2
PC	PC-15013	Nodule 3
PC	PC-15014	Nodule 4
PC	PC-15015	Nodule 5

9.8 Overview of the Applied Patient Root Query IOD

This section specifies in detail the applied attributes in the C-FIND Service Element of this supported SOP Class.

The Attribute Specific Character Set (0008,0005) is not included. The scanner will use predefined DICOM Character Sets depending upon the selected language. See table in section 7.

The values of key attributes are set in the search dialog where a list of key attributes is displayed.

Types of attributes used in Q/R Information Models (Column “Type”):

U: Unique Key Attribute

R: Required Key Attribute

O: Optional Key Attribute

Unique, required and options keys may be contained in the identifier of a C-FIND request.

Types of matching performed on key attributes (Column “Match”):

S: Single Value Matching

L: List of UID Matching

U: Universal Matching

W: Wild Card Matching

R: Range Matching

Q: Sequence Matching

Table 9.8-1: Key Attributes for Patient Root Q/R Information Model - C-FIND

Description	Tag	Type	Match	Note
Patient Level				
Patient’s Name	(0010,0010)	R	W, U	
Patient ID	(0010,0020)	U	W, U	
Patient’s Birth Date	(0010,0030)	O	R, S	
Patient’s Sex	(0010,0040)	O	U	
Number of Patient Related Studies	(0020,1200)	O	U	

Study Level				
Patient ID	(0010,0020)	U	S	
Study Instance UID	(0020,000D)	U	U	
Study Date	(0008,0020)	R	R, S	
Study Time	(0008,0030)	R	U	
Accession Number	(0008,0050)	R	S, U	
Study ID	(0020,0010)	R	S, U	
Referring Physician's Name	(0008,0090)	O	W, U	
Study Description	(0008,1030)	O	U	
Number of Study Related Series	(0020,1206)	O	U	
Number of Study Related Instances	(0020,1208)	O	U	
Series Level				
Patient ID	(0010,0020)	U	S	
Study Instance UID	(0020,000D)	U	S	
Series Instance UID	(0020,000E)	U	U	
Modality	(0008,0060)	R	S, U	
Series Number	(0020,0011)	R	U	
Composite Object Instance Level				
Patient ID	(0010,0020)	U	S	
Study Instance UID	(0020,000D)	U	S	
Series Instance UID	(0020,000E)	U	S	
SOP Instance UID	(0008,0018)	U	U	
Instance Number	(0020,0013)	R	U	

9.9 Overview of the Applied Patient Root/Study Root Retrieve IOD

This section specifies in detail the applied attributes in the C-MOVE Service Element of this supported SOP Class.

The Attribute Specific Character Set (0008,0005) is not included. The scanner will use predefined DICOM Character Sets depending upon the selected language. See table in section 7.

Unique Keys are contained in the identifier of C-MOVE. Required and Optional Keys are not contained in the identifier of C-MOVE.

If the level of retrieval is STUDY, the move destination is the local AE title, and the default listening port for C-STORE sub-operation is 104. If the level of retrieval is IMAGE, the move destination is local AE title suffixed with “_SCP”, and the default listening port for C-STORE sub-operation is 7878. When configuring AE titles and listening ports for move destination, attention should be taken to populate both of the above two entries.

Table 9.9-1: Key Attributes for Patient Root / Study Root Q/R Information Model - C-MOVE

Description	Tag	Note
Patient ID	(0010,0020)	Single value of Unique Key Attribute
Study Instance UID	(0020,000D)	Single value of Unique Key Attribute
Series Instance UID	(0020,000E)	Single value of Unique Key Attribute. Empty if retrieval level is STUDY
SOP Instance UID	(0008,0018)	List of UIDs. Empty if retrieval level is STUDY

Table 9.9-2: Accepted Presentation Contexts for C-STORE Sub-operation of C-MOVE

Abstract Syntax		Transfer Syntax	Role	Extended Negotiation
Name	UID			
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.830.10008.1.2.2 1.2.840.10008.1.2.4.57 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.90 1.2.840.10008.1.2.4.91 1.2.840.10008.1.2.4.92 1.2.840.10008.1.2.4.93	SCP	None

Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.830.10008.1.2.2 1.2.840.10008.1.2.4.57 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.90 1.2.840.10008.1.2.4.91 1.2.840.10008.1.2.4.92 1.2.840.10008.1.2.4.93	SCP	None
Comprehensive Structured Report Storage	1.2.840.10008.5.1.4.1.1.88.33	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.830.10008.1.2.2	SCP	None

Carestream Health, Inc.

150 Verona Street

Rochester, NY 14608

United States

© Carestream Health, Inc., 2015

Printed in the United States.

CARESTREAM is a trademark of Carestream Health.

PN 9J9407

