

# Wireless Tablet Upgrade.

SMART IMAGING SOLUTIONS

### Take Mobile Imaging to a New Level – with X-Factor Thinking.

With the CARESTREAM DRX-Transportable System/Universal Mobile Retrofit Kit, you can upgrade portable analog X-ray systems from almost any manufacturer to the power of DR wireless technology. Easy to install, this innovative kit helps you maximize your investment in your current mobile units.

The Universal Mobile Retrofit Kit is a non-invasive, flexible retrofit package. Self-powered and rechargeable, it works with a DR detector to transfer images instantaneously, so they can quickly be transferred to physicians for patient diagnosis and treatment planning. And, since you won't have to contend with tethers or cables during examinations, this is a system you'll use with ease and convenience in any healthcare facility.

### Advantages At a Glance

- Quickly and cost-effectively convert existing analog mobile systems to DR.
- It's easy to use and easy to upgrade.
- Wireless capability improves productivity and workflow.
- The system is self-powered and rechargeable.
- Choose from a wide range of available detector options.





# DRX-Transportable System/Universal Mobile

## The Power of DR Imaging, Easier and More Portable Than Ever.

- Upgrade existing equipment to DR.
- Wireless capability lets you go anywhere.
- Universal mounting hardware for the tablet PC is included.
- Choose from a variety of CARESTREAM DR Detector options:
  - DRX Plus** This versatile detector features the X-Factor and is available in multiple sizes, including a 25 x 30 small-format detector that's ideal for pediatric imaging.
  - Lux 35** Our uniquely designed glass-free detector includes ergonomic features for easier positioning and more comfortable patient exams.
  - Focus 35C and 43C** These provide economical ways to upgrade to full wireless DR technology.



## Product Specifications

### Key Components:

- Tablet PC: HP Elite X2 G4 Tablet
  - 8th Generation Intel® Core™ i5 Processor
  - 13.0 inch Multi-Touch Display with 3000 x 2000 resolution
  - 16GB of Memory
  - 512GB Solid-State Drive for Storage
  - 1.4 kg (3.1 lbs)
- DR Detectors (sold separately)



## Trusted Partner: CARESTREAM Smart Services

At Carestream, we'll be there when you need us. Our medical imaging solutions are backed by our worldwide support team of product and industry specialists – experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important – your patients.

## X-Factor Thinking

The X-Factor makes the most of your capital spending with configurable, modular components powered by the portable DRX Plus and Lux detectors. Use the detector on your retrofitted mobile units for your morning rounds and later, take it to a busy DRX-Evolution room for peak patient volumes. When the X-Factor lets you scale your solution to fit your workflow and budget, the possibilities are endless.



## Eclipse: Harness the Power of Decades of Innovation.

Carestream's Eclipse is the engine behind our innovative imaging software. It uses AI technology and proprietary algorithms to significantly increase the value of the entire imaging chain, from capture to diagnosis.

## Industry-Leading Performance: Meet Our Family of DR Detectors.

Our portfolio of DR detectors includes wireless, shareable, glass-free and fixed models to meet your demands based on exam type, detector size, dosage level and budget. Discover the right detector for your needs at [carestream.com/detectors](http://carestream.com/detectors).

[carestream.com](http://carestream.com)



"Rx only"

© Carestream Health, Inc., 2022. CARESTREAM is a trademark of Carestream Health. CAT 2000 003 03/22

