

Carestream

DRX-Excel System

XRS

# TOTAL FLEXIBILITY.

Fluoroscopy and Radiology  
in a Single System.



# A TIMELY TWO-IN-ONE SOLUTION.

The popularity of fluoroscopy for diagnostic imaging is surging – and with good reason. Fluoroscopy provides excellent image quality and supports a wide range of exams.

Now this technology has been optimized specifically for hospitals and imaging centers, with the CARESTREAM DRX-Excel System. It combines both fluoroscopy and general radiology capabilities in one compact unit. It also delivers accelerated workflow, high-resolution images and easy, convenient operation. *That's total flexibility.*

# AFFORDABLE VERSATILITY FOR HOSPITALS AND IMAGING CENTERS.

Meet the ideal blend of performance and economy. The DRX-Excel System gives you capabilities for both general radiography and fluoroscopy exams in a single system – at a very attractive price. And the system's versatility doesn't end there: it's also available in a choice of two models to match your needs and workflow:

- A conventional configuration, utilizing film or CR cassettes or our wireless DRX-Detector for radiology, and an image intensifier for the fluoroscopy exam.
- Or a DRF configuration, with a single, flat-panel digital detector for both fluoroscopy and radiology.

## Simply a Better Solution.

If you're shopping for a combination Fluoroscopy/RAD system, be sure to shop carefully. Compared to the DRX-Excel, various competing models have significant disadvantages. These include:

- Limited source-to-image distance limit of 150cm.
- Shorter tabletops.
- The need to manually rotate the tube.
- Lower table-weight limits.
- The absence of four-way table motion.
- Inability to adjust any lower than 65cm.

The DRX-Excel can transcend these limitations – providing an outstanding experience for both technologists and patients.




This easy-to-use console features a large color touch-screen and intuitive controls. It integrates all table movements, generator settings and APR functions into a single, efficient interface. Auto-positioning facilitates faster exam set-up times.

The optional DRX-1 System with an integrated console offers the power of the X-Factor and allows for all radiology exams to be done whether in table, cross-table or weight-bearing.



# FULLY-FEATURED FLEXIBILITY.

The DRX-Excel's feature set delivers imaging performance of the highest order. Ergonomic design and attention to detail help ensure patient comfort and user productivity. A focus on efficient system operation accelerates workflow. And high-resolution image capture and processing support excellent image quality and fast, accurate diagnoses.



## Keep a Close Watch on Your Patient.

A built-in video camera provides you with a live feed, so you can monitor the patient's position and overall status during the entire exam. The console's microphone lets you speak with your patient for instructions or reassurance whenever needed.



## Music To Your Patients' Ears.

The DRX-Excel (DRX-Excel Plus) can play musical selections to relax and entertain patients during their exam. You can select your own customized menu of music tracks during system installation. Plus, you can choose to have these music selections accompanied by soothing ambient lighting effects.

**Right for Today. Ready for Tomorrow.**

Carestream is ready to help you plan the most effective route to your X-ray imaging future. Our scalable equipment design and modular components mean high performance today, along with easy, affordable upgrades for years to come.

Let's plan out your equipment migration together – you'll gain the confidence that your current technology investment will continue to pay dividends well into your future.



## DYNAMIC TABLE CAPABILITIES.

The DRX-Excel offers an advanced, remote-controlled table with exceptional flexibility for a wide range of applications. Its ergonomics and efficiency optimize the entire examination process – with convenience for the operator and a low-stress, comfortable experience for the patient.

## HIGH-PERFORMANCE OPTIONS.

**Four-way table movement** facilitates easier patient positioning.

**Image-stitching station** combines multiple images into one for easier reading of spine and lower-limb exams.



Table tilts +90°/-30° to accommodate a wide range of exam types.



Table features a 200kg weight capacity to accommodate heavier patients.

[carestream.com/drx-excel](http://carestream.com/drx-excel)

"Rx only"

