



## Mobile X-ray Imaging is Really on the Move.

Now, X-Factor Thinking lets you take the DRX Detector to the field with this rugged, all-in-one solution. It includes the wireless detector, electronics, wireless access point and a tablet PC – all securely packed in a sturdy, portable carrying case. The DRX-Transportable System/Field Portable is ideal for applications such as the military, veterinary medicine or travel between nursing homes and homecare patients. You'll have all the performance of a DR imaging room – wherever you need it.



## ADVANTAGES AT A GLANCE

- Compact, protective carrying case for the DRX-Transportable makes it easy to transport to remote locations.
- Quickly and cost-effectively convert existing analog mobile systems to DR.
- The system is easy to use, with minimal training required.
- Wireless capability improves productivity and workflow.
- It's available with the DRX Plus 3543, DRX Plus 3543C and DRX 2530C Detectors.

# CARESTREAM DRX-Transportable System/Field Portable

## The Power of DR Imaging, More Portable Than Ever Before.

- Capture and view images immediately at the point of patient care.
- Lightweight, wireless DRX Detectors are easy to use in confined spaces.
- The system is self-powered and rechargeable.
- The rugged carrying case is water- and dust-resistant and easy to transport.



## Your Trusted Partner: CARESTREAM Smart Services

At Carestream, we'll be there when you need us. Our medical-imaging solutions are backed by our worldwide support team of product and industry specialists – experts who work closely with you to keep your equipment up and running at peak performance. By delivering a comprehensive scope of service offerings, we provide complete peace of mind over the life of your equipment, so you can focus on what's most important – your patients.

## Product Specifications

### Key Components:

- Thin, 454 x 368 x 48 mm (17.9 x 14.5 x 1.9 in.) electronics box containing the wireless access point, Ethernet switch, generator interface board and power supply
- Same type of battery used for electronics box and DRX detector; battery charges while plugged in
- Tablet PC
- Hardware to mount tablet PC to system
- Carrying case measuring 625 x 500 x 297 mm (24.6 x 19.7 x 11.7 in.)
- DRX Detectors\* – DRX Plus 3543, DRX Plus 3543C and DRX 2530C Detectors

## Right for Today. Ready for Tomorrow.

Carestream is ready to help you plan the most effective route to your X-ray imaging future. Our scalable equipment design and modular components mean high performance today, along with easy, affordable upgrades for years to come.

Let's plan your equipment migration together – you'll gain the confidence that your current technology investment will continue to pay dividends well into your future.

## A Community of Service and Support

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.

\*Detectors sold separately.



[carestream.com](http://carestream.com)

