

IHE Integration Statement		Date:	December 2009
Vendor	Product Name	Version	
Carestream Health	Image Suite	1.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles implemented	Actors Implemented	Options Implemented	
Portable Data for Imaging (PDI)	Portable Media Creator	None	
Internet address for Carestream Health IHE information: http://www.carestreamhealth.com/ihe.html			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.carestreamhealth.com/dicom.html		
Links to general information on IHE			
North America: http://www.rsna.org/IHE	Europe: http://www.ihe-europe.org	Japan: http://www.jira-net.or.jp/ihe-j	