

IHE Integration Statement		Date:	15 Aug 2003
Vendor	Product Name	Version	
Eastman Kodak Company	<i>Kodak DryView 8200 Laser Imager with Kodak PACS Link Systems</i>	5.2	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server	None	
Basic Security (SEC)	Secure Node	None	
Internet address for Kodak IHE information: www.kodak.com/go/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	www.kodak.com/go/dicom		
Links to general information on IHE			
North America: www.rsna.org/IHE	Europe: www.ihe-europe.org	Japan: www.jira-net.or.jp/ihe-j	