


Ultra-speed



## The intraoral film that delivers high-contrast and fine detail.

Ultra-speed dental film is a D-speed film that generates the high-quality, detailed images you need to make accurate diagnoses. Ideal for both manual and automated processing methods, Ultra-speed film is a versatile and easy-to-use solution for your practice.

### Enhanced image quality.

Because it accommodates such a wide range of processing and handling conditions, Ultra-speed film is a simple solution for your intraoral imaging needs. Exceptionally tolerant of process variations, this innovative film is designed to render excellent image quality.

### The standard in safety.

Ultra-speed film is also available with a ClinAsept barrier feature to help reduce the risk of infection. The built-in barrier protects the film packet from contaminants, while

its innovative design frees up your time and ensures peace of mind. Simply open the barrier envelope, place the film into the cup and you're ready to start processing.

Available in multiple sizes and formats to meet all of your intraoral imaging needs:

- Available in sizes 0-4.
- Select from flexible, easy-to-open Super Poly-Soft vinyl packets or sturdier paper packets for size 2.
- Bite-wing packets – paper packets with bite-wing tabs attached.
- ClinAsept barrier packets for added safety.
- For optimum results, process your radiographs using Carestream Dental chemicals.

## Ordering Information for Ultra-speed

Size	Description	Code	Packets/Package	REF No
Size 0	Super Poly-Soft packets	DF-54	100	122 8840
Size 0	Super Poly-Soft packets (double film)	DF-53	100	122 8931
Size 0	Super Poly-Soft packets with ClinAsept barrier	DF-54C	75	144 5360
Size 1	Paper packets	DF-56	100	127 3747
Size 1	Paper packets (double film)	DF-55	100	127 3721
Size 2	Super Poly-Soft packets	DF-58	150	165 8194
Size 2	Super Poly-Soft packets (double film)	DF-57	130	175 3664
Size 2	Paper packets	DF-58	150	149 1737
Size 2	Paper packets (double film)	DF-57	150	149 1752
Size 2	Bite-wing paper packets	DF-40	50	839 3043
Size 2	Super Poly-Soft packets with ClinAsept barrier	DF-58C	100	171 7131
Size 2	Super Poly-Soft packets with ClinAsept barrier (double film)	DF-57C	100	898 7877
Size 3	Bite-wing paper packets	DF-42	100	129 6771
Size 4	Occlusal paper packets	DF-50	25	166 6163

For more information, call **800.933.8031** or visit [carestream.com](http://carestream.com).