

IHE Integration Statement		Date:	October, 2009
Vendor	Product Name	Version	
Carestream Health	CARESTREAM RIS	11.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	DSS/Order Filler PPS Manager	Appointment Notification	
Patient Information Reconciliation (PIR)	DSS/Order Filler PPS Manager	None	
Charge Posting (CHG)	DSS/Order Filler	None	
Audit Trail and Node Authentication (ATNA)	Secure Node DSS/Order Filler PPS Manager	None	
Consistent Time (CT)	Time Client	None	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-referencing Consumer	None	
Retrieve Information for Display (RID)	Display Information Source	none Summary of All Reports Summary of Radiology Reports	
Internet address for Carestream Health IHE information: ihe.carestreamhealth.com			
Links to Standards Conformance Statements for the Implementation			
HL7	www.hl7.org		
DICOM	dicom.carestreamhealth.com		
Links to general information on IHE			
North America: www.rsna.org/IHE	Europe: www.ihe-europe.net	Japan: www.jira-net.or.jp/ihe-j	