


READYMATIC / READYPRO



Carestream Dental film and chemicals: the perfect combination

READYMATIC and READYPRO dental chemicals have been specially developed for use with INSIGHT and Ultra-speed dental films. The result: consistently high resolution images.



When our dental chemicals are used to process our dental film, you'll get:

- Outstanding image quality for accurate diagnosis
- Superior clarity and resolution
- Compatibility with all leading dental film processors
- Simplified workflow – ready-to-use formulas mean no starter solution, mixing or dilution

Superb Results with Minimal Hassle

READYMATIC dental chemicals

READYMATIC dental chemicals are designed for use in all automatic roller-transport dental x-ray film processors. They deliver clear, consistent results for all intraoral and extraoral radiographs. Plus, they:

- Provide excellent film transportability and sensitometry
- Are specifically formulated to limit residue on rollers and in tanks
- Will keep for up to two years, unopened and stored between 5–30 °C (40–85 °F)
- Are available in convenient chem packs or as separate solutions

READYPRO Twin Pack

READYPRO developer and fixer are specially formulated to keep rollerless processors running. Plus, they:

- Come packaged in pre-measured bottles; all you have to do is open and pour
- Will keep for up to two years, unopened and stored between 5–30 °C (40–85 °F)
- Are available in handy twin packs

Other Carestream Dental Chemicals

Ask us about our wide variety of other chemistry products including X-OMAT chemicals for large automatic processors, Rapid Access dental chemicals for quick access to images and GBX chemicals for manual processing.

Ordering Information

Product Name	Description	REF No
READYMATIC Chem Pack	2x5 L	860 6873
READYMATIC Developer and Replenisher	2x5 L	528 5929
READYMATIC Fixer and Replenisher	2x5 L	528 5937
READYPRO Twin Pack	4x2 L	860 6881

For more information, call **800.933.8031** or visit carestream.com.