

IHE Integration Statement		Date:	Oct 2004
Vendor	Product Name	Version	
Eastman Kodak Company	<i>Kodak DirectView</i> Web Software	5.2	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display	None	
Consistent Presentation of Images (CPI)	Image Display Print Composer	None	
Access to Radiology Information (ARI)	Image Display	None	
Basic Security (SEC)	Secure Node Image Display	None	
Internet address for Kodak IHE information: www.kodak.com/go/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	www.kodak.com/go/dicom		
Links to general information on IHE			
North America: www.rsna.org/IHE	Europe: www.ihe-europe.org	Japan: www.jira-net.or.jp/ihe-j	