



# **IHE Integration Statement CARESTREAM Vue PACS v12.1**

**Part # 9J8036  
2015-06-29**

<b>Integrating the Healthcare Enterprise (IHE) Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Carestream Health, Inc	CARESTREAM Vue Archive (Vue for Clinical Collaboration Platform), CARESTREAM Vue PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry.	CARESTREAM Vue PACS, Vue Archive v12.1.0.0	2015-06-29
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed in this table.			
<b>Note:</b> This document is intended for use by Carestream customers.			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None	
Imaging Object Change Management (IOCM)	Image Manager/Image Archive	None	
Access to Radiology Information (ARI)	Image Manager/Image Archive Image Display Report Repository External Report Repository Access	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None	
Cross-enterprise Document Reliable Interchange (XDR)	Document Recipient Document Source	None	

<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Cross-enterprise Document Sharing (XDS.b)	Document Consumer Document Repository Document Registry Document Source	Multiple document
Cross-enterprise Document Sharing for Imaging (XDS.b-1)	Imaging Document Consumer Imaging Document Source	Set of DICOM Instances, Multiple documents None
Patient Demographic Query HL7 V3 (PDQv3)	Patient Demographics Consumer	None
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None
Evidence Documents (ED)	Image Display Image Manager/Image Archive	None
Cardiology Evidence Documents (ED-CARD)	Image Display Image Manager/Image Archive	None
Simple Image and Numeric Report (SINR)	Report Repository External Report Repository Access	None
Mammography Image (MAMMO)	Image Manager/Image Archive Image Display	None
Key Image Note (KIN)	Evidence Creator Image Display Image Manager/Image Archive	None

<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None
Audit Trail and Node Authentication (ATNA)	Secured Application	None
Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web Content None None None
<b>Link to Carestream Health IHE information:</b> <a href="http://ihe.carestream.com">ihe.carestream.com</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://dicom.carestreamhealth.com">dicom.carestreamhealth.com</a>	
<b>Links to General Information on IHE</b>		
In North America: <a href="http://www.rsna.org/ihe.aspx">www.rsna.org/ihe.aspx</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>



Carestream Health France  
1, rue Galilée  
93192 NOISY-LE-GRAND CEDEX  
FRANCE



**Carestream**



Carestream Health, Inc.  
150 Verona Street  
Rochester, NY 14608  
United States

© Carestream Health, Inc. 2015.

Created in the USA.

CARESTREAM is a trademark of Carestream Health, Inc.

Pub No. 9J8036  
Rev A